

## Steris/Amsco Harmony LED 585 Surgical Lights - Refurbished



Stock Photo – Options May Vary

The Steris Harmony 585 LED surgical lighting system is designed to provide high intensity illumination in almost any surgical setting. The 24" light head is easily maneuvered and long-life LEDs boast an expected life of 30,000 hours.

The Harmony LEDs produce up to 14,800 ft candles (160k lux) of cool, pure white, natural light in a pattern that can be adjusted from 7" to 11". The wall mounted control center is included and allows for intensity adjustment at the control center or on the light heads themselves.

The Steris Harmony LED surgical light can come in a single or dual head configuration.

The Steris Harmony 585 surgical light system come fully refurbished and include a comprehensive warranty.

- Cutting Edge LED Technology
- Pure White, Natural Light
- Excellent Depth of Illumination (26")
- Color temperature 4400k
- 360° Spring/Swivel Arm Support
- Intensity: 14,800 ft candles (160,000 lux)
- Wall Mounted Intensity Controller
- Color Rendering Index: 97
- Focusable Light Head
- Excellent Shadow Control
- Long Life LED Bulbs
- Completely Refurbished

### Specifications:

- Intensity Pattern (large): 11"
- Intensity Pattern (small): 7"
- Clearance (rotation): 16' (approx)
- LED Life: Approx 30,000 hrs
- Intensity: Up to 14,800 Ft Candles
- CRI: >97
- Color Temp: 4400 K
- Adjustable Light Pattern: 7" – 11"
- Diameter of Light Head: 24"
- Power Requirement (volts/Hz):
  - 100-240V, 50/60 Hz
- Power Requirements (amps):
  - 6A (dual head)
  - 3A (single head)

*\*\*Brochure and images are a general overview of the product.  
Actual product options may vary.*