

## Stryker Rugged Ambulance Stretcher - Refurbished



*Stock Photo – Options May Vary*

Stryker Rugged series Ambulance Stretcher is a perfect choice for emergency transportation. The lightweight, yet durable construction makes the stretcher easy to handle and navigate. The one handed operation of the backrest and side rails simplifies use. The basic model has a 500lbs weight capacity, but other models in the series can handle up to 650 lbs.

The break-away head section helps in tight areas, reducing the stretcher length by to 62”, making the Stryker Rugged series one of the most maneuverable stretchers on the market today.

The Stryker Rugged stretchers come completely refurbished with a comprehensive warranty and a new pad.

- Transport at load height capability
- Positive action height adjustment
- Easy-to-use release handle design
- Color-coded controls
- Lightweight rugged aluminum construction
- 500 lbs Weight Capacity
- One-hand release, breakaway head section
- One-hand release, fold down side rails
- One-hand release, infinite positioning, pneumatically assisted backrest
- Completely Refurbished

### Specifications:

- Overall Length: 62” - 80.5”
- Overall Width: 23”
- Weight: 81 lbs
- Patient Surface Width: 26”
- Height Range: 13.5” – 37.5”
- Weight Capacity: 500 lbs
- Trend/Rev Trend: ± 18°
- Side Rails: 14" x 54"
- Wheel Diameter: 6”

\*\*Frame types may vary

*\*\*Brochure and images are a general overview of the product. Actual product options may vary.*