



[www.venturemedical.com](http://www.venturemedical.com)

Venture Medical ReQuip, Inc.  
6008 Bonacker Drive  
Tampa, FL 33610  
800-627-3215 | 813-740-8640  
Fax: 800-466-1251  
[sales@venturemedical.com](mailto:sales@venturemedical.com)

## ValleyLab Force 2 Electrosurgical Unit - Refurbished



- Digital Display
- Multipurpose ESU
- High Power Output
- Accurate and Consistent Cutting
- Three Blend Modes
- REM Fault Indicator
- Low Voltage Coag Mode
- Monopolar Foot-Switch Included
- Cart (Optional)
- Completely Refurbished

The Valleylab Force 2 Electrosurgical Unit (ESU) is a powerful multipurpose ESU that is suitable for a wide range of procedures.

Features include: 3 Preset Blend Modes and Enhanced Safety with the REM Return Electrode Monitoring System. The REM system is vital in helping to prevent electrode burns.

This unit also features 300 watts of max cut power, 120 watts in coag, and 70 watts in bipolar. The easy to read digital display shows the settings and audible alarms alert the user to problems.

The ValleyLab Force 2 comes completely refurbished with a comprehensive warranty and includes a Monopolar Foot-switch, Powercord and Manual.

### Specifications:

- Monopolar Cut: 300 Watts
- Monopolar Coag: 120 Watts
- Bipolar: 70 Watts
- Blend 1: 250 Watts
- Blend 2: 200 Watts
- Blend 3: 150 Watts
- Low Voltage Coag: 99 Watts
- Height (unit only): 8"
- Width (unit only): 13"
- Length (unit only): 21"
- Weight (unit only): 23 lbs
- Power Requirements: 120V, 6A Max

# VENTUREMEDICAL.COM

© 2007-2018 Venture Medical ReQuip, Inc.