

ValleyLab Force FX Electrosurgical Unit - Refurbished



Stock Photo – Options May Vary

The ValleyLab Force FX Electrosurgical Unit (ESU) provides a smoother cut through all tissue types and automatically adjusts to changes in tissue, while maintaining the needed power for the surgeon.

Features include: 3 Cutting Modes (Low Cut, Pure Cut, and Blend for Cutting w/ Hemostasis), Instant Response Technology, and Power Efficiency Rating (PER) of 98 for Accurate and Consistent Cut Performance.

The ValleyLab Force FX comes completely refurbished with a comprehensive warranty and includes a Monopolar Foot-switch, Powercord and Manual.

- Instant Response Technology
- Auto Adjust to Tissue Changes
- Accurate and Consistent Cutting
- 3 Cutting Modes
 - Low Cut
 - Pure Cut
 - Blend
- Power Efficiency Rating (PER): 98
- Accurate and Consistent Cut Performance
- 4 Coag Modes
- Monopolar Foot-Switch Included
- Cart (optional)
- Completely Refurbished

Specifications:

- Monopolar Cut: 300 Watts
- Monopolar Coag: 120 Watts
- Bipolar: 70 Watts
- Height (unit only): 11.1 cm (4-3/8")
- Width (unit only): 35.6 cm (14")
- Length (unit only): 43.9 cm (17")
- Weight (unit only): 8.1 kg (18 lbs)
- Power Requirements: 120V, 7A Max

***Brochure and images are a general overview of the product. Actual product options may vary.*