

## Wy'East Totalift II Stretcher Chair - Refurbished



Stock Photo – Options May Vary

\*\* Please Note: Upper Side Rails and Contoured Headrest are Optional

The Wy'East Totalift II Stretcher Chair is an extremely versatile stretcher that is perfect for a wide variety of departments within your medical facility.

The Totalift Stretcher features a pneumatic assist fowler (backrest), hydraulic height adjustment, central brake and steering system, side rails, O2 tank well, and sliding patient surface for easy transfers.

The Wy'East Totalift II comes completely refurbished with a comprehensive warranty and a new pad.

- Transfer Slide Top
- Gas Spring Assist Backrest
- ± Height, + Trendelenburg
- Tuck-Away Side Rails
- 4 Wheel Locking
- Built-in O2 Tank Holder
- 400 lbs Weight Capacity
- Includes New 2" Pad
- Completely Refurbished

### Specifications:

- Overall Length: 71½"
- Overall Width: 27"
- Height Range: 23" - 36"
- Weight Capacity: 400 lbs
- Fowler Angle: 0° - 80°
- Trendelenburg: 20°
- Wheel Diameter: 4"
- Pad Thickness: 2"

*\*\*Brochure and images are a general overview of the product.  
Actual product options may vary.*