

Zoll M-Series Defibrillator Biphasic | Pacing | AED - Refurbished



Stock Photo – Options May Vary

The Zoll M-Series defibrillator is one of the smallest and most lightweight defibrillators available in its class. This unit features biphasic waveform, non-invasive pacing and an AED advisory mode. Operating the system is simple and user-friendly with straightforward controls and easy to use menus.

Included with the unit: Multi-Function Cable, ECG Cable and Leadwires, Set of Adult Therapy Pads, Printer Paper and New Battery.

Optional accessories and functions include: 12 Lead ECG, SpO2, NIBP, ETCO2 (Mainstream or sidestream), Carrying Case and External Paddles.

The Zoll M-Series comes completely refurbished with a comprehensive warranty.

- Biphasic Waveform Technology
- Pacing Function
- AED Function
- Synchronized Cardioversion
- 3 Lead EKG
- Onboard Thermal Printer
- Hands Free Multi-Function Cable
- Built-In Battery Charger
- Small, Lightweight Design
- Carry Handle
- Many Optional Functions Available
- Monitor Displays Selected and Delivered Energy
- Completely Refurbished

Specifications:

- Max Joules: 360j
- Biphasic up to 200j
- Charge Time: Less than 7 seconds
- Display Size (Screen): 14.2 cm (5.6")
- Battery: Rechargeable, Sealed Lead Acid
- Recorder Paper: 80mm, Thermal
- Height (unit only): 17.3 cm (6.8")
- Width (unit only): 26.2 cm (10.3")
- Depth (unit only): 20.8 cm (8.2")
- Weight (unit only): ~5.23 kg (~11.5 lbs)
- Power Requirement: 120v / 2 amps

***Brochure and images are a general overview of the product. Actual product options may vary.*