


BOVIE® DERM 102 HIGH FREQUENCY DESICCATOR



The Bovie® DERM 102 High Frequency Desiccator offers ten watts of monopolar power, adjustable in 1/10th watt increments, to desiccate and fulgurate tissue. The power can be adjusted using the 3-Button Power Control Handpiece A902 (Patent #8,998,899), or using the rotary dial on the unit. The unit features a large, bright, LED display for excellent visibility across the procedure room. If foot control activation is preferred while using the DERM 102, just plug in the A803 Bovie Footswitch (sold separately) and activate. No separate handpiece is required.

For the practitioner who wants to also use bi-polar electrosurgery along with high frequency desiccation, the DERM 102 High Frequency Desiccator provides 10 watts of bipolar power, adjustable in 1/10th watt increments. To utilize this feature, the A827V* bipolar cable, a bipolar forcep* and footswitch activation are required.

The DERM 102 is packaged ready to use. The unit comes with a total of 18 disposable dermal tips. By the time the practitioners utilizes all 18 electrodes, they will have paid for the unit (twice), and generated income for their practice.

BOVIE® DERM 102

HIGH FREQUENCY DESICCATOR

PRODUCTS INCLUDED

 <p>DERM 102</p>	
 <p>A902</p>	 <p>A804, A805, A806, A807, A806DE, A807DE</p>
 <p>A837</p>	 <p>A910</p>

Style # DERM 102

▼ Products Included

- DERM 102** (1) high frequency desiccator
A902 (1) handpiece
A910 (2) non-sterile, disposable handpiece sheaths
A804 (5) non-sterile single-use sharp dermal tips
A805 (2) sterile single-use sharp dermal tips
A806 (5) non-sterile single-use blunt dermal tips
A807 (2) sterile single-use blunt dermal tips
A806DE (2) non-sterile DERM-Elite™ Premium blunt dermal tips
A807DE (2) sterile DERM-Elite Premium blunt dermal tips
A837 (1) wall mount kit
 110 VAC Hospital-grade power cord*
 User's guide CD (**MC-55-229-001**)*

Dimensions

Height: 208mm (8.2")
 Width: 184mm (7.2")
 Depth: 69mm (2.7")
 Weight: < 1.4kg (< 3 lbs.)

Output Power

Up to 10 watts (0-10 watts in 1/10th watt increments on both monopolar and bipolar)

Line Voltage

100-240 VAC

Line Frequency: 50 – 60 Hz

Main Current: 320 mA Max

Warranty: FOUR YEARS

* Not Pictured