

Drager Narkomed GS Anesthesia Machine - Refurbished



Stock Photo – Options May Vary

The Drager Narkomed GS Workstation offers an excellent combination of quality and performance, friendly design and advanced features. The Narkomed GS is a perfect answer for medical facilities that want to upgrade to a more modern machine, yet remain budget conscious.

The GS features: AV2+ Electronic Ventilator, Absorber, PEEP, and Volume. Pressure and O₂ Monitors are built in. Additional features include Fresh Gas System, Scavenger System, and O₂ Flush.

The Narkomed GS comes completely refurbished with a comprehensive warranty.

- Budget Conscious
- Reliable
- AV2+ Electronic Ventilator
- Absorber
- PEEP
- Volume, Pressure and O₂ Monitors Built In
- Fresh Gas System
- Scavenger System
- O₂ Flush
- 2 or 3 Vaporizer Manifold
- 3 Drawers for Storage
- Alarms
- Battery Backup System
- Completely Refurbished

Specifications:

- Overall Dims: 36" W x 53.25"H x 27" D
- Weight: 375 lbs
- Tidal Volume Monitoring:
 - Range: 0.2-50l
 - Resolution: 0.1 L
 - Accuracy: 10% or 0.1 L, whichever is greater
- Airway Pressure Monitoring:
 - Range: -10-125 cm H₂O
 - Resolution: 1cm H₂O
 - Accuracy: ±3cm H₂O or ±10 of reading, whichever is greater
- Oxygen Monitoring:
 - Range: 10-100 Vol%
 - Resolution: 1 Vol%
 - Accuracy: ±3 Vol%

***Brochure and images are a general overview of the product.
Actual product options may vary.*