



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Tuttnauer 1730M Valueklave Autoclave - Refurbished



- Manual Autoclave
- Economical
- Overheat Protection
- Automatic Shutoff
- Drying Cycle
- Includes Trays
- 7" Diameter Chamber
- Front Facing Drain Valve
- Completely Refurbished
- Comprehensive Warranty

The Tuttnauer 1730M Valueklave is one of the best selling tabletop automatic/electronic autoclaves in the industry. The 1730M features a 7" diameter chamber and is both reliable and economical. It is the perfect choice for your basic sterilization needs.

The Tuttnauer 1730M features a fast cycle time of only 11 minutes from a hot start (standard, unwrapped) and offers an efficient drying cycle for wrapped instruments. Its compact size also is perfect for tight spaces and the standard electrical plug means it can be operated from a normal electrical outlet. The unit is also equipped with 3 trays.

The Tuttnauer 1730M Valueklave comes completely refurbished with a comprehensive warranty and includes trays.

Specifications:

- Chamber Dims: 7" Dia x 13" Depth
- Overall Depth: 17.9"
- Overall Width: 17.4"
- Overall Height: 12"
- Number of Trays: 3
- Trays: 11.6" x 4.7" x 0.8" (DxWxH)
- Volume: 2 Gallons
- Electrical Requirements:
 - 120 Volt, 11 Amps
- NRTL Certification Listed

VENTUREMEDICAL.COM

© 2007-2025 Venture Medical ReQuip, Inc.