



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Tuttnauer 3870M Autoclave - Refurbished



The Tuttnauer 3870M is one of the largest volume tabletop autoclaves in the industry. The 3870M features a 15" diameter chamber and is both reliable and economical. It is the perfect choice for larger sterilization needs without upgrading to a hospital-grade unit.

The Tuttnauer 3870M features a fast cycle time for its capacity of only 23 minutes from a hot start (standard, unwrapped) and offers an efficient drying cycle for wrapped instruments. The unit is also equipped with 2 trays.

The unit has an auto shut-off at the end of both the sterilization and drying cycles, and a double safety locking device prevents accidental openings of the door when pressurized. The drain valve is conveniently located on the front of the unit for easy draining.

The Tuttnauer 3870M comes completely refurbished with a comprehensive warranty and includes trays.

- Manual Autoclave
- Large Capacity Chamber (22 gal)
- Deep 30" Chamber
- 15" Diameter Chamber
- Economical
- Overheat Protection
- Automatic Shutoff
- Drying Cycle
- Includes Trays
- Front Facing Drain Valve
- Completely Refurbished
- Comprehensive Warranty

Specifications:

- Chamber Dims: 15" Dia x 30" Depth
- Overall Depth: 34.5"
- Overall Width: 26"
- Overall Height: 20.6"
- Number of Trays: 2
- Volume: 22 Gallons
- Electrical Requirements:
 - 220/230 Volt, 13.6 Amps
- NRTL Certification Listed