

Stryker SM104 Stretcher - Refurbished



Stock Photo – Options May Vary

The Stryker SM104 stretcher features mobility, simplicity, and versatility. Fully equipped and completely adaptable, the M-Series allows you to choose the mobility solution that best fits your needs. Supremely mobile and user-friendly, these stretchers greatly reduce the physical strain on clinicians without sacrificing patient care and comfort.

Additional features include: power-washable design with 700 lbs weight capacity, 4 integrated IV receptacles, litter-mounted O2 holder, and glide-away siderails.

The SM104 features all standard features plus fifth wheel steering: uses a retractable fifth wheel to provide superior maneuverability and cornering.

The Stryker SM104 stretcher comes completely refurbished with a comprehensive warranty and a new pad.

- Gas Spring Assist Backrest
- Glide-away Side Rails
- Central Brake and Steer
- \pm Height
- \pm Trendelenburg
- Built-in O2 Tank Holder
- 700 lbs Weight Capacity
- Includes New 4" Pad
- Completely Refurbished

Specifications:

- Overall Length: 84"
- Overall Width: 33.5"
- Patient Surface Width: 26" x 75.25"
- Weight capacity: 700 lbs
- Height range: 20.75" - 34.5"
- Knee Gatch: 40°
- Trend/Rev Trend: \pm 16°
- Side Rails: 14" x 57.75"
- Caster Diameter: 8"
- Pad Thickness: 4"

***Brochure and images are a general overview of the product.
Actual product options may vary.*