

Tuttnauer 1730M Valueklave Autoclave - Refurbished



Stock Photo – Options May Vary

- Manual Autoclave
- Economical
- Overheat Protection
- Automatic Shutoff
- Drying Cycle
- Includes Trays
- 7" Diameter Chamber
- Front Facing Drain Valve
- Completely Refurbished
- Comprehensive Warranty

The Tuttnauer 1730M Valueklave is one of the best-selling tabletop automatic/electronic autoclaves in the industry. The 1730M features a 7" diameter chamber and is both reliable and economical. It is the perfect choice for your basic sterilization needs.

The Tuttnauer 1730M features a fast cycle time of only 11 minutes from a hot start (standard, unwrapped) and offers an efficient drying cycle for wrapped instruments. Its compact size also is perfect for tight spaces and the standard electrical plug means it can be operated from a normal electrical outlet. The unit is also equipped with 3 trays.

The Tuttnauer 1730M Valueklave comes completely refurbished with a comprehensive warranty and includes trays.

Specifications:

- Chamber Dims: 7" Dia x 13" Depth
- Overall Depth: 17.9"
- Overall Width: 17.4"
- Overall Height: 12"
- Number of Trays: 3
- Trays: 11.6" x 4.7" x 0.8" (DxWxH)
- Volume: 2 Gallons
- Electrical Requirements:
 - 120 Volt, 11 Amps
- NRTL Certification Listed

***Brochure and images are a general overview of the product.
Actual product options may vary.*