

Tuttnauer EZ9 Fully Automatic Autoclave - Refurbished



Stock Photo – Options May Vary

- Automatic Autoclave
- Closed door active drying system
- Automatic Shutoff
- Includes Trays
- 9” Diameter Chamber
- Front Facing Drain Valve
- Completely Refurbished
- Comprehensive Warranty

The EZ9 with the Simplicity of One Touch Design! With the simplicity of one touch design all your sterilization and drying needs are fulfilled. The Tuttnauer EZ9 fully automatic autoclave will fill, sterilize, exhaust and dry at the touch of a button. Closed door active drying system to maintain sterility and ensure efficient drying of packs and pouches.

Additional Tuttnauer EZ9 Features: Hepa filter 0.2 µm air filter (replaceable). All program parameters can be customized and stored, Automatic shut off at the end of both the sterilization and dry cycles, and Audible and visual cycle interruption alert. International certifications and comes with a 2-year parts and labor warranty.

Specifications:

- Chamber Dims: 9” Dia x 18” Depth
- Overall Depth: 21.5”
- Overall Width: 20”
- Overall Height: 14.4”
- Number of Trays: 2 Full & 2 Half
- Trays: 16.3” x 6.7” x 0.8” (DxWxH)
- Volume: 5 Gallons
- Electrical Requirements:
 - 120 Volt

***Brochure and images are a general overview of the product.
Actual product options may vary.*