

## Luvis M200 Mobile LED Surgical Light - New



*Stock Photo – Options May Vary*

The Luvis M200 LED Surgical Light is the perfect light for nearly any procedure. The M200 features exceptional quality with energy efficient operation.

The 54 LED'S in each light head produce up to 120,000 lux in a pattern that can be adjusted from 7.8" to 11.8". The Luvis M200 Surgical Light also features a control panel that is synchronized a sub-handle for the user to utilize the control panel easily during surgery.

The Luvis M200 surgical light can come in a single, dual or mobile configuration and includes a 2-year mfg. warranty.

- Excellent Illumination
- Excellent Dilution Effect
- Ergonomic & Friendly Design
- Detachable & Autoclavable Handle
- No UV Emission & Cool Light
- Excellent Depth of Illumination
- 4,300k Color Temperature
- Intuitive and Easy Way of Control
- Rotating Central Knob
- Intensity: 120,000 lux
- Color Rendering Index: 95
- Focusable Light Head
- Lightweight Suspension System
- Various Installation Options
  - Ceiling, Mobile, Dual Ceiling Mount

### Specifications:

- LEDs: 54
- CRI: 95
- Color Temperature: 4,300k
- Head Size: 19.6"
- Focus Size: 7.8" (Circle)/11.8" (Circle)
- Focal Distance: 1M
- Intensity: Max. 120,000lux
- Input Power: 100-240Vac, 60Hz
- Power Consumption: 85W
- LED Life Span: 50,000hr

*\*\*Brochure and images are a general overview of the product. Actual product options may vary.*