

## Skytron 6002 Elite Operating Table - Refurbished



Stock Photo – Options May Vary

The Skytron 6002 Elite operating table is a modern multi-purpose surgical table and a great choice for almost any surgical setting. The 6002 features include wired hand control, 180° tabletop rotation, hydro-electric brakes, C-Arm compatibility, and a 600 lbs. weight capacity in normal position.

The Elite operating table has all standard functions ± back, ± leg, ± height, ± tilt, ± flex, trend/rev trend, Return to Level, and comes with North American side rails which fit most accessories.

The Skytron 6002 Elite comes completely refurbished with a comprehensive warranty and a new 2' pad set.

- C-Arm Compatible Table Top
- Wired Hand Control
- 600 lbs Weight Capacity
- 180° Tabletop Rotation
- 90° Backrest Elevation
- Return to Level
- Electro-Hydraulic Floor Locks
- Stainless Steel Base and Column
- North American Side Rails
- Swivel Casters
- Includes New 2" Pad Set

### Specifications:

- Overall Length: 75"
- Overall Width: 22"
- Patient Surface Width: 20"
- Height Range: 28" - 43"
- Weight Capacity: 600 lbs
- Trend/Rev Trend: ± 30°
- Back Rest: + 90° / - 40°
- Leg Rest: 105°
- Lateral Tilt: 30°
- Flex and Reflex Capable
- Power: 120 Volts
- Pad Thickness: 2"

*\*\*Brochure and images are a general overview of the product.  
Actual product options may vary.*