



HCPCS  
E0600  
A7000  
A7002



18610  
18615

Two configurations  
(AC or AC/DC) using  
one efficient design



## VacuMax™ Plus and VacuMax™ Go Suction Machines

| Item # | Product Description                  | Product HCPCS       | UOM  |
|--------|--------------------------------------|---------------------|------|
| 18610  | VacuMax™ Plus Suction Machine        | E0600, A7000, A7002 | 1/ea |
| 18615  | VacuMax™ Go Portable Suction Machine | E0600, A7000, A7002 | 1/ea |

- Easy-to-turn adjustable flow regulator knob conveniently placed directly below gauge for easy use and control on both models.
- Both models include a 800cc suction canister, 6' suction tube, 10" suction tubing, hydrophobic filter, elbow connector and manual.
- The Drive 18615 VacuMax™ Go Portable Suction Machine has all the same features and benefits that the 18610 VacuMax™ Plus provides, as well as an internal rechargeable, replaceable battery that lasts up to 60 minutes per charge.
- The 18610 may be converted into the 18615 with the purchase and provider installation of the 18615-BATTERY.
- The 18615 may be charged while in use and requires less than 8 hours to reach a full charge.
- The Drive 18615 VacuMax™ Go is ideal for providing suction in the home or on the go, which comes complete with rechargeable battery, carry bag and a DC power cord for the car.
- The Drive VacuMax™ vacuum adjustments allow up to 150-530 mmHg and a free flow up to 25 lpm.

| accessories   |  | UOM   |
|---------------|--|-------|
| 18615-ACORD   | AC Power Cord for the 18610 and 18615        | 1/ea  |
| 18615-DCORD   | DC Power Cord for the 18615                  | 1/ea  |
| 18615-BATTERY | Replaceable VacuMax™ Go Rechargeable Battery | 1/ea  |
| 18615-BAG     | Carry Bag for VacuMax™ Suction Machine       | 1/ea  |
| 18600-FILTER  | Hydrophobic Filter for Suction Machines      | 3/bx  |
| 22330         | 800cc Disposable Suction Canister Kit        | 1/ea  |
| 22330-12      | 800cc Disposable Suction Canister Kit        | 12/cs |
| 18600-KITN    | Universal Suction Tubing and Filter Kit      | 1/bx  |
| 610-12B       | 800cc Disposable Suction Canister            | 12/cs |
| 610-48BP      | 800cc Disposable Suction Canister            | 48/cs |

### specifications

|                             |                |
|-----------------------------|----------------|
| Warranty on Suction Machine | Limited 2 Year |
| Warranty on Battery         | 6 months       |