



## Advantages

- EMR Connectivity**  
Scale has the ability to interface with a computer, EMR software or other electronic device via various Health o meter® Professional Connectivity Solutions
- Calculates Body Mass Index (BMI)**  
Scale functions include the ability to calculate the patient's Body Mass Index (BMI)
- Value-Added Features**  
Operating functions (LB/KG Conversion, LB/KG Lock Out, BMI, Zero Out/Tare, Hold/Release, Reweigh, Recall, Auto Zero, Auto Off) included at no additional cost
- Hospital Grade Construction**  
Durable steel structure with large sturdy stretcher enables patients to be weighed accurately and comfortably

# Health o meter®

P r o f e s s i o n a l



## Specifications

- Capacity: 400 lb / 181 kg
- Resolution: 0.2 lb / 0.1 kg
- Connectivity via USB
- 3/4" LCD display
- Canvas Size: 70½" (w) x 30" (d) / 1791 mm (w) / 762 mm (d)
- Product Footprint: 70½" (w) x 46½" (d) x 60½" / 1791 mm (w) x 1181 mm (d) x 1537 (h)
- Product Weight: 114lb / 52kg
- Functions: LB / KG Conversion, LB / KG lock Out, Body Mass Index (BMI), Zero Out / Tare, Hold / Release, Reweigh, Auto Zero, Auto Off
- (4) Wheels for ease of movement, 1 locking wheel for safety
- Power Source: 100-240V adapter included, 6 AA batteries (not included)
- 2 Year Limited Warranty