



PresSion™ 4-Chamber Garment Measurement Chart



**651-4 PresSion Pump,
Part no 4307**

Lower Extremity		Part No.	Length (in.)	Circumference (in.**)		
Size	Type		Inseam	Ankle	Knee	Thigh/Groin
Small	Full Leg	43197	23	21	25	30
Small Wide	Full Leg	43198	26	24	29	35
Small Extra Wide	Full Leg	43199	24	28	35	40
Medium	Full Leg	43200	29	22	27	31
Medium Wide	Full Leg	43201	29	25	31	35
Medium Extra Wide	Full Leg	43202	23	29	36	40
Large	Full Leg	43203	34	21	27	31
Large Wide	Full Leg	43204	36	22	31	40
Half	Half Leg	43205	16.5	22	28	–
Half Wide	Half Leg	43206	19	24	30	–
Upper Extremity		Part No.	Length (in.*)	Circumference (in.**)		
Size	Type		Axilla to Fingertips	Axilla	Elbow	Wrist
Medium	Arm	43207	26.5	22	19	13
Large	Arm	43208	32	23	20	14

Lower Extremity		Part No.	Length (cm.)	Circumference (cm.**)		
Size	Type		Inseam	Ankle	Knee	Thigh/Groin
Small	Full Leg	43197	58.4	53.3	63.5	76.2
Small Wide	Full Leg	43198	66	61	73.7	88.9
Small Extra Wide	Full Leg	43199	61	71.1	88.9	101.6
Medium	Full Leg	43200	73.7	55.9	68.6	78.8
Medium Wide	Full Leg	43201	73.7	63.5	78.7	88.9
Medium Extra Wide	Full Leg	43202	71.1	73.7	91.4	101.6
Large	Full Leg	43203	86.4	53.3	68.6	78.7
Large Wide	Full Leg	43204	91.4	55.9	78.7	101.6
Half	Half Leg	43205	41.9	55.9	71.1	–
Half Wide	Half Leg	43206	48.3	61	76.2	–
Upper Extremity		Part No.	Length (cm.*)	Circumference (cm.**)		
Size	Type		Axilla to Fingertips	Axilla	Elbow	Wrist
Medium	Arm	43207	67.3	55.9	48.3	33
Large	Arm	43208	81.3	58.4	50.8	35.6

*Length corresponds to the entire sleeve length, from axilla to fingertips.

**Indicates the maximum circumference measurement of each body part that will fit into the corresponding sleeve size shown.



PresSion™ 8-Chamber Garment Measurement Chart



**652-8 PresSion Pump,
Part no 4309**

Lower Extremity		Part No.	Length (in.)	Circumference (in.**)		
Size	Type		Inseam	Ankle	Knee	Thigh/Groin
Small	Full Leg	43209	23	21	25	30
Small Wide	Full Leg	43210	26	24	29	35
Small Extra Wide	Full Leg	43211	24	28	35	40
Medium	Full Leg	43212	29	22	27	31
Medium Wide	Full Leg	43213	29	25	31	35
Medium Extra Wide	Full Leg	43214	23	29	36	40
Large	Full Leg	43215	34	21	27	31
Large Wide	Full Leg	43216	36	22	31	40
Upper Extremity		Part No.	Length (in.*)	Circumference (in.**)		
Size	Type		Axilla to Fingertips	Axilla	Elbow	Wrist
Medium	Arm	43217	26.5	22	19	13
Large	Arm	43218	32	23	20	14
Medium	Arm & Shoulder	43219	26.5	22	19	13
Large	Arm & Shoulder	43220	32	23	20	14

Lower Extremity		Part No.	Length (cm.)	Circumference (cm.**)		
Size	Type		Inseam	Ankle	Knee	Thigh/Groin
Small	Full Leg	43209	58.4	53.3	63.5	76.2
Small Wide	Full Leg	43210	66	61	73.7	88.9
Small Extra Wide	Full Leg	43211	61	71.1	88.9	101.6
Medium	Full Leg	43212	73.7	55.9	68.6	78.8
Medium Wide	Full Leg	43213	73.7	63.5	78.7	88.9
Medium Extra Wide	Full Leg	43214	71.1	73.7	91.4	101.6
Large	Full Leg	43215	86.4	53.3	68.6	78.7
Large Wide	Full Leg	43216	91.4	55.9	78.7	101.6
Upper Extremity		Part No.	Length (cm.*)	Circumference (cm.**)		
Size	Type		Axilla to Fingertips	Axilla	Elbow	Wrist
Medium	Arm	43217	67.3	55.9	48.3	33
Large	Arm	43218	81.3	58.4	50.8	35.6
Medium	Arm & Shoulder	43219	67.3	55.9	48.3	33
Large	Arm & Shoulder	43220	81.3	58.4	50.8	35.6

*Length corresponds to the entire sleeve length, from axilla to fingertips.

**Indicates the maximum circumference measurement of each body part that will fit into the corresponding sleeve size shown.



DJO, LLC | A DJO Global Company
T 800.592.7329 **D** 760.727.1280 **F** 800.936.6569
 1430 Decision Street | Vista, CA 92081-8553 | U.S.A.
DJOglobal.com/chattanooga

Together in Motion™