

DB-3000

Surgical Table Accessories



The Rigid Ergonomic Drainage Bag System is designed to fit any North American surgical table and is compatible with all suction systems. Efficient fluid containment for hysteroscopy and cystoscopy procedures.

Specifications

Weight	20lb
Warranty	2 Year

Features:

- Integral tissue retention and filter with flexible hose.
- Non-Flexible Support Frame
- Integral rail clamps and one bag included.
- Sterile & Non-sterile bags available.
- Attaches level w/the surgical table side rails.

Materials:

- Non-Flexible, solid stainless steel support frame
- Drainage Bag
- Integral Rail Clamps

