

# GP-2800

## Gel Body Positioners

The Contoured Ankle Pads are made from viscoelastic gel that offers the highest quality and most comfortable positioners available with superior anti-shear protection. The fluid-resistant gel helps protect against skin tissue breakdown and other types of costly post-operative damage. Kova Gel™ Pads are Latex Free, Fluid Resistant, Reusable and Covered with a Membrane-Like Coating for Easy Clean Up. Perfect for decreasing shear movement and reducing post-operative complications while providing maximum protection and comfort for the patient.



### Features:

- 4" x 7.5" x 3" - 1 pair
- Anti-shear protection
- Provides support, stability, and elevation
- Latex free and fluid resistant
- MRI Compatible
- Fluid Proof
- Reusable - Easy to Clean
- Disinfectable

### Materials:

- Kova Gel viscoelastic polyurethane gel
- Vinyl membrane covering

### Specifications

Weight	7lb
Pieces	1 pair
Warranty	1 Year