

GP-3300

Gel Body Positioners

The Candy Cane Style Legholder Post Pads are made from viscoelastic gel that offers the highest quality and most comfortable positioners available with superior anti-shear protection. The fluid-resistant gel helps protect against skin tissue breakdown and other types of costly post-operative damage. They are Latex free, fluid resistant, reusable and covered with a membrane-like coating for easy clean up.

Features:

- Anti-shear protection
- 2 leg-holder post gel pads
- Provides support, stability, and elevation
- Latex free and fluid resistant
- MRI Compatible
- Fluid Proof
- Reusable - Easy to Clean
- Disinfectable

Materials:

- Kova Gel viscoelastic polyurethane gel
- Vinyl membrane covering



Specifications

| | |
|----------|---------|
| Weight | 5lb |
| Pieces | 1 piece |
| Warranty | 1 Year |