

# GP-3350

## Gel Body Positioners

The Candy Cane Style Legholder Ankle Straps are made from viscoelastic gel that offers the highest quality and most comfortable positioners available with superior anti-shear protection. The fluid-resistant gel helps protect against skin tissue breakdown and other types of costly post-operative damage. They are Latex free, fluid resistant, reusable and covered with a membrane-like coating for easy clean up. Patient safety straps are made specifically for the medical industry to help and insure the patients stability and proper positioning.

### Features:

- Anti-shear protection
- 12"L x 1.5"W x 0.375" Thick - 2pc
- Provides support, stability, and elevation
- Latex free and fluid resistant
- MRI Compatible
- Fluid Proof
- Reusable - Easy to Clean
- Disinfectable

### Materials:

- Kova Gel viscoelastic polyurethane gel
- Vinyl membrane covering



### Specifications

Weight	7lb
Pieces	1 piece
Warranty	1 Year