

GP-3580

Gel Pillow Positioners



The Multi-Positioner Pillow Gel Filled Sand Bag is made from viscoelastic gel that offers the highest quality and most comfortable positioners available with superior anti-shear protection. The fluid-resistant gel helps protect against skin tissue breakdown and other types of costly post-operative damage. Kova Gel™ Pads are Latex Free, Fluid Resistant, Reusable and Covered with a Membrane-Like Coating for Easy Clean Up. Perfect for decreasing shear movement and reducing post-operative complications while providing maximum protection and comfort for the patient.



Features:

- Anti-shear protection
- 5" x 14" x 3"
- Provides support, stability, and elevation
- Latex free and fluid resistant
- MRI Compatible
- Fluid Proof
- Reusable - Easy to Clean
- Disinfectable

Materials:

- Kova Gel viscoelastic polyurethane gel
- Vinyl membrane covering

| Specifications | |
|----------------|---------|
| Weight | 6lb |
| Pieces | 1 piece |
| Warranty | 1 Year |