

GP-3750

Gel Pillow Positioners

The Multi-Positioner Pillow Gel Filled Sand Bag is made from viscoelastic gel that offers the highest quality and most comfortable positioners available with superior anti-shear protection. The fluid-resistant gel helps protect against skin tissue breakdown and other types of costly post-operative damage. Kova Gel™ Pads are Latex Free, Fluid Resistant, Reusable and Covered with a Membrane-Like Coating for Easy Clean Up. Perfect for decreasing shear movement and reducing post-operative complications while providing maximum protection and comfort for the patient.



Features:

- Anti-shear protection
- 5" x 20" x 2¾"
- Provides support, stability, and elevation
- Latex free and fluid resistant
- MRI Compatible
- Fluid Proof
- Reusable - Easy to Clean
- Disinfectable

Materials:

- Kova Gel viscoelastic polyurethane gel
- Vinyl membrane covering

Specifications	
Weight	6lb
Pieces	1 piece
Warranty	1 Year