

GP-5350

Gel Body Positioners

The Patient Restraint Strap with Kova Gel™ Pad dimensions are 3½" x 20" x ¾" and 2" x 114" Strap. The Kova Gel™ Patient Safety Straps secures to the side of the table. The Strap is made of nylon and the pad is made of viscoelastic polyurethane gel, this helps to reduce tissue trauma. They are Latex free, fluid resistant, reusable and covered with a membrane-like coating for easy clean up. Patient safety straps are made specifically for the medical industry to help and insure the patients stability and proper positioning. They are used in the surgery room on Surgical Tables, Armboards, or Beds, and also used on stretchers for a safer transfer.

Specifications

Weight	7lb
Pieces	1 piece
Warranty	1 Year

Features:

- Anti-shear protection
- Gel Pad: ½" x 20" x ¾"
- Strap: 2" x 114"
- Latex free and fluid resistant
- MRI Compatible
- Fluid Proof
- Reusable - Easy to Clean
- Disinfectable

Materials:

- Kova Gel viscoelastic polyurethane gel
- Vinyl membrane covering
- Nylon Strap

