

DERM 101

**Bovie®**



# MAKING THE OFFICE PRACTITIONER MORE PRODUCTIVE



HIGH-FREQUENCY DESICCATOR

## ECONOMICAL

Affordably priced to minimize the facility's expense and maximize their ROI.

## AMPLE POWER

Ten (10) watts of monopolar power which are adjustable in 1/10<sup>th</sup> watt increments.

## SIMPLICITY

The thoughtfully designed unit allows Practitioners the ability to perform minor skin procedures with an easy-to-operate device, to increase their practice revenue.

## EQUIPPED

Packaged complete and ready to use. By the time the practitioners utilize all 18 electrodes they will have paid for the unit, twice.

## RELIABILITY

The Derm 101 is backed by Bovie's standard four (4) year warranty.



DERM 101  
INCLUDES

			
DERM 101 Unit	A902 3-Button Handpiece	A804 Sharp Dermal Tips, Non-Sterile (5ea)	A805 Sharp Dermal Tips, Sterile (2ea)
			
A806 Blunt Dermal Tips, Non-Sterile (5ea)	A807 Blunt Dermal Tips, Sterile (2ea)	A806DE Derm-Elite™ Premium Blunt Dermal Tip, Non-Sterile (2ea)	A807DE Derm-Elite™ Premium Blunt Dermal Tip, Sterile (2ea)
			
A910 Disposable Handpiece Sheath, Non Sterile (2ea)	A837 Wall Mount Kit (1ea)	Medical Grade Power Cord	User Manual (CD Format)

DERM 102  
& ACCESSORIES

The Bovie® Derm 102 has an additional Bipolar mode offering up to 10 watts of power in 1/10th watt increments.





		
A803 Footswitch*	Reusable Bipolar Forceps* †	A827V Bipolar Forcep Cable*

\*Sold Separately † More styles available

Additional accessories available—contact your Bovie representative for details

DERM 101 and DERM 102 Technical Specifications...

<b>Dimensions and Weight:</b>	<b>Mains Connection:</b>
Height 2.7" (69mm) Width 7.2" (184mm) Length 8.2" (208mm) Weight < 1.4 kg (< 3 lbs)	Main Voltage: 100 – 240 VAC ± 10% Main Frequency: 50 – 60 Hertz Main Current: 320mA Max. Power Consumption: 28 VA
<b>Safety:</b>	 
Basic Construction: In accordance with EN 60601-1 Mode of Operation: Intermittent operation Output Type: TYPE BF Duty Cycle: 10sec on / 30sec off Main Fuses: T 400mAL, 250V	

TECHNICAL  
SPECIFICATIONS



Bovie Medical Corporation • 5115 Ulmerton Road • Clearwater, FL 33760 • www.boviemed.com • sales@boviemed.com  
U.S. Phone 1-800-537-2790 • Fax 1-800-323-1640 • Int'l Phone +1-727-384-2323 • Fax +1-727-347-9144

Bovie® is a registered trademark of Bovie Medical Corporation