



[www.venturemedical.com](http://www.venturemedical.com)

Venture Medical ReQuip, Inc.  
6008 Bonacker Drive  
Tampa, FL 33610  
800-627-3215 | 813-740-8640  
Fax: 800-466-1251  
[sales@venturemedical.com](mailto:sales@venturemedical.com)

## Conmed 2450 Electrosurgical Unit - Refurbished



- Digital Display
- Bipolar and Monopolar
- Multipurpose ESU
- High Power Output
- Accurate and Consistent Cutting
- Programmable
- ARM Grounding Fault Indicator
- Accurate and Consistent Cut Performance
- Dynamic Response Technology
- Monopolar Foot-Switch Included
- Optional Bipolar Foot Switch and Cart
- Completely Refurbished

The Conmed 2450 Electrosurgery generator (ESU) is a state-of-the-art unit with an easy to use operator interface. Smooth cutting, precise bipolar, and ARM technology make the Conmed 2450 safe and effective, with less tissue damage.

The Conmed System 2450 Electrosurgical Unit features 2 Cutting Modes: Pure Cut for scalpel cuts and Blend for hemostasis while cutting. Also, the unit has coagulation for both dissection and hemostasis, and precise Bipolar desiccation with minimal lateral tissue damage. The Bipolar setting includes both Micro and Macro modes.

Dynamic Response Technology reacts to tissue changes to reduce tissue damage, and the on-board Automatic Return Monitor (ARM) detects faults in the patient grounding circuit to help prevent pad site burns.

Included with the unit: Monopolar Foot Switch and Users Manual. Options include Bipolar Foot Switch and Cart.

The Conmed 2450 ESU comes completely refurbished with a comprehensive warranty.

### Specifications:

- Monopolar Cut (Max): 300 Watts
- Monopolar Coag: 120 Watts
- Bipolar: 70 Watts
- Blend Cut: 200 Watts
- Hi Blend Cut: 200 Watts
- Standard Coag: 120 Watts
- Spray Coag: 80 Watts
- Micro Coag (BiPolar): 70 Watts
- Macro Coag (Bipolar): 70 Watts
- Power Requirements: 120V
- Weight: 32 Lbs (unit only)
- Dims: 14.4"W x 21.5"D x 7"H