

## Conmed 5000 Electrosurgical Unit - Refurbished



Stock Photo – Options May Vary

The Conmed 5000 Electrosurgery generator (ESU) combines technology and versatility providing the ability to fine-tune the ESU to specific procedures and surgeon's preferences.

Unique modes on the System 5000™ include the Fluids Mode, for an additional boost of energy in procedures taking place in a fluid environment, and Lap Mode, which limits open circuit voltages to below 2700 volts in all modes and reduces the potential risk of inadvertent burns from capacitive coupling.

The Automatic Return Monitor (ARM™) is the most sensitive pad monitoring system available. It combines automatic feedback circuitry and visual displays to provide the most sensitive contact quality monitoring system available today. The ARM™ function continuously monitors the safety of the patient to ground System pad interface with a dual foil ground pad. It disarms the ESU and alerts the staff if a pad fault condition occurs.

Included with the unit: Monopolar Foot Switch and Users Manual. Options include Bipolar Foot Switch and Cart.

The Conmed 5000 ESU comes completely refurbished with a comprehensive warranty.

- Digital Display
- Bipolar and Monopolar
- Multipurpose ESU
- High Power Output
- Accurate and Consistent Cutting
- 9 Programmable Settings
- ARM Grounding Fault Indicator
- Accurate and Consistent Cut Performance
- 3 Advanced Specialty Modes
- Monopolar Foot-Switch Included
- Optional Bipolar Foot Switch and Cart
- Completely Refurbished

### Specifications:

- Monopolar Cut (Max): 300 Watts
- Monopolar Coag: 120 Watts
- Bipolar: 50 Watts
- Blend 1: 200 Watts
- Blend 2: 200 Watts
- Blend 3: 200 Watts
- Spray Coag: 80 Watts
- Micro Coag (BiPolar): 50 Watts
- Macro Coag (Bipolar): 90 Watts
- Power Requirements: 115V, 5.3A

*\*\*Brochure and images are a general overview of the product. Actual product options may vary.*