



www.venturemedical.com

Venture Medical ReQuip, Inc.  
6008 Bonacker Drive  
Tampa, FL 33610  
800-627-3215 | 813-740-8640  
Fax: 800-466-1251  
sales@venturemedical.com

## Luvis M300 LED Dual Surgical Light - New



- Excellent Illumination
- Ergonomic & User-Friendly Design
- Detachable & Autoclavable Main Handle
- No UV Emission & Cool Light
- Excellent Depth of Illumination
- 4,300k Color Temperature
- Intuitive and Easy to Use Controls
- Rotating Central Knob
- Intensity: 160,000 lux
- Color Rendering Index: 95
- Focusable Light Head
- Lightweight Suspension System
- Various Installation Options
  - Ceiling, Mobile, Dual Ceiling Mount

The Luvis M300 LED Surgical Light is the perfect light for any procedure. The M300 features exceptional quality with energy efficient operation.

The 84 LED'S in each light head produce up to 160,000 lux in a pattern that can be adjusted from 7.8" to 11.8". The Luvis Surgical Light also features a control panel that is synchronized with sub-handle for the user to approach the control panel easily during surgery.

The Luvis M300 surgical light can come in a single, dual or mobile configuration and includes a 2 year mfg. warranty.

### Specifications:

- LEDs: 84
- CRI: 95
- Color Temperature: 4,300k
- Head Size: 23.6"
- Focus Size: 7.8" (Circle)/11.8" (Circle)
- Focal Distance: 1M
- Intensity: Max. 160,000lux
- Input Power: 100-264Vac, 60Hz
- Power Consumption: 125W
- LED Life Span: 50,000hr

**VENTUREMEDICAL.COM**

© 2007-2025 Venture Medical ReQuip, Inc.