



Eschmann **operating table**
accessories



Quality & performance
you can trust, completely



ESCHMANN

Nov 2014



Eschmann **operating table** accessories


In terms of quality and reliability Eschmann operating tables are acknowledged as being second to none, the very best in medical engineering technology for over 100 years – our assurance for your complete confidence.

The Eschmann philosophy is one of working closely with the people that matter the most - the surgical clinicians and operating theatre professionals that use our world class products. This close cooperation allows us to develop a range of operating table accessories that consistently exceed the exacting demands required in today's operating theatre environment.

This brochure showcases the current range of Eschmann operating table accessories, allowing surgeons and medical professionals to undertake almost every type of procedure.




Exceptional support and performance



Specialist traction accessories



Easy anatomical positioning



Secure specialist positioning

Accessories detailed in this brochure can be used on all operating tables unless otherwise stated. Suggested clamp requirements are stated against each accessory. Please ensure that the correct clamps are ordered to suit the table for which the accessories are intended.



General Surgery

Clamps, arm, foot, bariatric and general support accessories

4



Maxillofacial and Ophthalmics

Range of modular patient positioning accessories

12



Bariatric Surgery

Range of accessories capable of supporting high patient weights

13



Gynaecology and Urology

Range of gynaecology and urology positioning equipment

14



Orthopaedic and Trauma

Range of specialist and modular orthopaedic accessories

16



Spinal and Neurosurgery

Specialist skull clamp and head support systems

23



Patient Positioning Aids

Range of modular patient positioning accessories

27



Ancillary Products

Replacement handsets and cables for T Series tables

33



General Accessories

Accessory fixing clamps

Rotary clamps will accept all standard 16mm diameter accessory posts. Some of the clamps shown are only used on certain accessories and are detailed here for ease of identification.

Please ensure the correct clamp is ordered to suit the table for which the accessories are intended.



Blade clamp
Attaches at table rail ends
UK/Euro sidebar 81-273-49
US Sidebar 81-273-52

Easy-lock blade clamp
UK sidebar fitting 81-269-78
Euro sidebar fitting TA-020-1187
US sidebar fitting 81-269-81


Drop handle clamp
UK/Euro sidebar fitting TA-020-1084
UK/US sidebar fitting TA-020-1085

Rotary anti drift clamp
for standard 16mm posts
UK sidebar fitting TA-020-1082
Euro sidebar fitting TA-020-1086
US sidebar fitting TA-020-1087

Rotary accessory clamp
for posts up to 18mm
UK sidebar 81-270-45
Euro sidebar 81-270-52
US sidebar 81-270-48



Radial clamp for 16-20mm posts - UK	R10-307-UK
Radial clamp for 16-20mm posts - Euro	R10-307
Radial clamp for 16-20mm posts - Universal	R10-308
Radial clamp - 'hook on' - UK	R10-305-UK
Radial clamp - 'hook on' - Euro (not shown)	R10-305
Radial clamp - 'hook on' - US (not shown)	R10-305-US



Accessory clamp for 20mm posts - UK	R10-301-UK
Accessory clamp for 20mm posts - Euro	R10-301



Flexible anaesthetic screen 81-273-04
Requires 2 drop handle clamps (TA-020-1084/5)
or 2 rotary anti drift clamps (TA-020-1082/6/7)



Anaesthetic carbon fibre screen with extensions and integral UK clamp R10-313-UK
(Height range: 530mm - 860mm
Full width: 1700mm)
As above with integral Euro clamp R10-313



Narrow arm table with trigger, integral clamp and pad (660mm x 140mm)
UK - Complete with 50mm Pad 81-270-29
Euro - Complete with 50mm Pad 81-270-27
UK - Complete with 80mm Pad 81-270-39
Euro - Complete with 80mm Pad 81-270-42
Spare 50mm Pad 81-271-06
Spare 80mm Pad 81-271-08



Anaesthetic screen with sleeve 81-286-18
(550mm x 480mm)
Requires 1 drop handle clamp (TA-020-1084/5)
or 1 rotary anti drift clamp (TA-020-1082/6/7)
Spare sleeve for anaesthetic screen 81-288-12



Narrow arm table with contoured pad 81-300-19
(630mm x 210mm)
Spare contoured pad (580mm x 210mm) 81-302-13
Standard flat pad (580mm x 210mm) TA-050-1030
Requires or 1 drop handle clamp (TA-020-1084/5)
or 1 rotary anti drift clamp (TA-020-1082/6/7)



Narrow arm table with integral clamp and pad (660mm x 140mm)
UK - Complete with 50mm Pad 81-270-25
Euro - Complete with 50mm Pad 81-270-19
UK - Complete with 80mm Pad 81-270-32
Euro - Complete with 80mm Pad 81-270-37
Spare 50mm Pad 81-271-06
Spare 80mm pad 81-271-08



Eschmann **operating**table accessories



Light, hand and arm surgery table **81-270-24**
(800mm x 360mm, with 40mm mattress)

Requires 2 easy-lock blade clamps **81-269-78** (UK)
TA-020-1187 (Euro)



Armboard with ball joint release and integral clamp and pad **R10-387-UK-H**
(610mm x 180mm)- UK

Aa above with integral clamp - Euro **R10-387-H**
Spare pad for armboard **R10-382-H**



UK - Rectangular arm table **81-820-20**
Euro - Rectangular arm table **81-820-30**
Denyer - Rectangular arm table **81-820-40**
US - Rectangular arm table **81-820-50**

X-ray translucent with pad, incorporating telescopic support and direct on bracket. (690mm x 490mm, maximum support height 1110mm)

Spare pad for arm table **81-822-13**
UK - Spare Direct on Bracket **81-820-43**
Euro - Spare Direct on Bracket **81-820-45**
Denyer - Spare Direct on Bracket **81-820-47**
US - Spare Direct on Bracket **81-820-53**




Carbon fibre hand and arm surgery table complete with 25mm mattress, integral UK clamps and supporting leg **R10-396-UK**
(815mm x 520mm)

As above with integral Euro clamp **R10-396**
Spare pad for hand and arm table **R10-397**



Elbow support with integral UK clamp **R10-366-UK**
Elbow support with integral US clamp **R10-366-US**
Elbow support with Integral Euro Clamp **R10-366**



Curved leg/arm supports **R10-552**
Pad for curved leg/arm supports **R10-553**



Shoulder rests 81-273-39

For use with patients up to 225kg
Complete with 2 blade clamps (81-273-49 UK/Euro)




Leg abduction support (vein boards) TA-060-1011

For use with T10, T20 and T30 tables
Also available (not shown) leg abduction support (81-378-62) for J and MR series



Wristlet 81-446-13

Requires 1 drop handle clamp
(TA-020-1084/5)



Perspex arm support 81-368-23

Right angled, one only (300mm x 181mm x 102mm)



Foot stocks (pair) 81-378-38

(210mm x 130mm x 120mm, with 75mm wide recess)



Step stands - Set of 2 86-250-18

(458mm x 300mm x 115mm)



Pair of safety rails with integral UK clamp R10-531-UK

(620mm x 290mm)

Pair of safety rails with integral Euro clamp R10-531

Cover for safety rails (pair) R10-531H



Adjustable cotsides - each (730mm x 310mm) 81-462-77

Suitable for UK, Euro and US sidebars



Eschmann **operating**table accessories



Tibia support - (275mm x 135mm diameter) **R10-367**

Suitable Clamps:
R10-301 (Euro rail) for 20mm posts
R10-301UK (UK rail) for 20mm posts
R10-308 (Universal rail) for 16-20mm posts



Knee stabiliser (Post length 330mm) **81-273-16**

Requires 1 blade clamp 81-273-49 (UK/Euro)
81-273-52 (US)

Spare pad for knee stabiliser **81-273-27**



Buttock support **81-412-15**

Tubular, with pad (300 x 120mm diameter)
Requires 1 drop handle clamp (TA-020-1084/5)

Spare pad for buttock support **81-416-14**



Pelvis Support **81-404-13**

Tubular with pad (300mm x 120mm diameter)
Requires 1 drop handle clamp (TA-020-1084/5)

Spare pad for pelvis support **81-416-14**



Laminectomy support with 50mm pad **81-569-80**

(550mm x 260mm)
Requires 2 rotary anti drift clamps
TA-020-1082/6/7

Spare pad (50mm) **81-569-99**



Laminectomy support pad (montreal mattress) **81-569-13**
(980mm x 495mm)



Adjustable/rotatable lateral support - narrow **R10-560**
(200mm x 80mm) with straight and contoured pads
Requires clamp R10-301/-UK

Spare Narrow support pad - contoured **R10-737**
Spare Narrow support pad - straight **R10-894**



Adjustable lateral support - wide **R10-565**
(200mm x 80mm) with straight and contoured pads
Requires clamp R10-301/-UK

Spare Wide support pad - contoured **R10-738**
Spare Wide support pad - straight **R10-895**



Divided Leg Section Plus
Divided leg section dimensions: 700mm x 289mm
Maximum patient weight capacity: 325kg

UK sidebar with 50mm Mattress	TA-030-3120
Euro sidebar with 50mm Mattress	TA-030-3122
US sidebar with 50mm Mattress	TA-030-3124
Denyer sidebar with 50mm Mattress	TA-030-3126

Universal Footplate - (complete with clamps)
Universal Footplate dimensions: 309mm x 309mm

To fit UK/Euro clamps	TA-030-3142
To fit US clamps	TA-030-3144
To fit Denyer clamps	TA-030-3145



Footrest/extension UK Sidebars TA-020-3037

Complete with 50mm mattress. For use with T Series tables only (545mm x 330mm)

Footrest/extension UK Sidebars TA-020-3038
Complete with 80mm mattress. For use with T Series tables only

Footrest/extension US Sidebars (bare) TA-020-3045

Footrest/extension Euro Sidebars (bare) TA-020-3004

Requires 2 drop handle clamps (TA-020-1084/5).

Spare pads:

50mm mattress for footrest/extension	T20-431-2094
80mm mattress for footrest/extension	T20-432-0001

Also available (not shown) footrest/extension with pad and accessory clamps (460mm x 240mm) for RX tables (81-376-33) and MR and J Series tables (81-377-14).
Spare pad (80mm) (81-386-21)

Also available


Carter Braine arm support with pad Overall height 400mm Requires 1 drop handle clamp TA-020-1084/5	81-290-10
Spare pad for Carter Braine (81-290-10)	81-292-15



Lateral brace kit 81-270-35

Kit comprises 1 circular and 2 rectangular supports
Requires 3 blade clamps 81-273-49 (UK/Euro) 81-273-52 (US)

Rectangular support with pad	81-270-75
Circular support with pad	81-270-70
Replacement circular pad	110584
Replacement rectangular pad	110585



X-ray cassette tray TA-080-1090

End loading with flexible handle for size of cassette up to 650mm x 380mm. Trigger handle (1.2m) for loading.
Suitable for use with T Series operating tables

Also available X-ray cassette tray suitable for MR, J3, RX tables 81-827-52 (not shown)



Lateral support, (150mm x 210mm) Adult 81-492-16

Complete with pad and clamp (UK/Euro)

Spare pad - Adult 81-494-10

Lateral support, (140mm x 160mm) Child 81-504-19

Complete with pad and clamp (UK/Euro)

Spare pad - Child 81-506-13



Eschmann **operating**table accessories



Angular UK accessory rail with clamp R10-359-UK
(Bar length 250mm, pole length 230mm)
Angular Euro accessory rail with clamp R10-359



UK extension rail with integral clamp R10-491-UK
Euro extension rail with integral clamp R10-491



Transfusion pole, with integral UK clamp 81-463-06
(Height range: 590mm - 1000mm)
Transfusion pole, with integral Euro clamp 81-463-17



Infusion stand with integral rotatable UK clamp R10-400-UK
(Height range: 985mm - 1550mm)
Infusion stand with integral rotatable Euro clamp R10-400




Angled transport handles with integral UK clamps R10-456-UK
Angled transport handles with integral Euro clamps R10-456
Transport handles with integral UK clamps R10-457-UK
Transport handles with integral Euro clamps R10-457



T20/30/50 foot controller TA-040-4003
Height and Trend movements only

Accessory trolleys



General accessory trolley

The general accessory trolley has been developed to carry a range of table accessories including the T Series shoulder positioner and articulating leg sections.

Trolley complete with accessory box (UK sidebar)	81-582-57
Trolley complete with accessory box (Euro sidebar)	81-582-68
Trolley complete with accessory box (Denyer sidebar)	81-582-79
Spare accessory box	81-582-82



T20/T30 shoulder positioner adaptor for General accessory trollies 81-582-57 and 81-582-68	81-582-94
T20/T30 shoulder positioner adaptor for T Series head and leg trolley (TA-040-1120)	81-582-71

Head and leg flap accessory trolley

T Series head and leg flap accessory trolley
TA-040-1120





Maxillofacial and Ophthalmics



Ophthalmic/ENT headrest with pad 81-466-40
Spare contoured pad TA-050-1006



Ophthalmic Head Section with square top (bare) 81-465-27
Spare rectangular pad - 50mm 81-092-57
Spare rectangular pad - 80mm 81-092-73



Head stabiliser 'Dixey Rubens' 81-375-28
(500mm x 270mm x 100mm)
150mm diameter head recess



Multi adjustable headrest R10-230-F
Fits directly into Eschmann T Series tables

Bariatric accessories

Divided Leg Section Plus and universal footplates

Divided Leg Section Plus for T Series - UK - complete with 50mm mattress	TA-030-3120
Universal Footplates with 50mm Mattress - UK	TA-030-3142



325kg

Maximum Patient weight



600mm x 105mm

Standard width extender - UK (50mm pad)	TA-020-2033
Standard width extender - UK (80mm PRS pad)	TA-020-2034
Standard width extender - Euro (50mm pad)	TA-020-2046
Standard width extender - Denyer (Pad not included)	TA-020-2062
For use with T20, T30, RX, MR and J series tables	

Spare 50mm pad TA-050-1096



350mm x 105mm

Short width extender - UK (50mm pad)	TA-020-2035
Short width extender - UK (80mm PRS pad)	TA-020-2036
Short width extender - Euro (50mm pad)	TA-020-2045
Short width extender - Denyer (Pad not included)	TA-020-2056

For use with T20, T30, RX, MR and J series tables

Spare 50mm pad TA-050-1097



363kg

Maximum Patient weight

Bariatric direct placement legholders	81-269-07
Requires easy-lock blade clamps 81-269-78 (UK) TA-020-1187 (Euro)	
Spare pads for bariatric direct placement legholders (pair)	81-269-84
Direct placement legholder cart	81-272-45



Bariatric thigh straps - integral UK clamps (Pad - 800mm x 180mm, strap length 1870mm)	81-272-84
Bariatric thigh straps - integral Euro clamps	81-272-87
Bariatric thigh straps - integral US clamps	81-272-90
Bariatric thigh straps - integral Denyer clamps	81-272-93

Bariatric surgery kit - 81-150-23

Comprises: width extenders TA-020-2033 & 2035, bariatric direct placement legholders 81-269-07, 2 x easy-lock blade clamps 81-269-78 and thigh straps 81-272-84.



Gynaecology and Urology

159 kg

Maximum
Patient weight



- | | |
|---|-----------|
| Direct placement legholders (includes UK/
Euro clamps) | 81-272-10 |
| Spare pads for Direct placement legholders
(pair) | 81-272-15 |
| Direct placement legholder cart | 81-272-45 |



- | | |
|--|-----------|
| Pair of articulating leg sections
Complete with pads and boots, only
compatible with T Series tables | 81-620-03 |
| Boot liner for articulating leg sections | 81-272-15 |
| Replacement upper pad for articulating legs | 81-620-25 |
| Replacement lower pad for articulating legs | 81-620-36 |



- | | |
|--|-----------|
| Paediatric Stirrup Package
(includes Junior Paediatric Stirrups 81-271-09 and Little
Paediatric Stirrups 81-271-24)
Requires 2 blade clamps (81-273-49 UK/Euro, 81-273-52 US) | 81-271-29 |
| Junior Paediatric Stirrups | 81-271-09 |
| Junior Paediatric Clamshell Pads (pair) | 81-271-21 |
| Little Paediatric Stirrups | 81-271-24 |
| Little Paediatric Clamshell Pads (pair) | 81-271-27 |



- | | |
|---|-----------|
| Lithotomy Supports with sleeves & straps
(Gel pads sold separately)
Requires 2 rotary anti drift clamps (TA-020-1082/6/7) | 81-250-15 |
| Pair of large gel straps for lithotomy supports | 81-256-27 |
| Pair of small gel straps for lithotomy supports | 81-256-35 |



- | | |
|--|-----------|
| Lloyd Davies lithotomy supports (pair)
with pads and straps
Requires 2 rotary anti drift clamps
(TA-020-1082/6/7) | 81-420-17 |
| Spare pads (pair) | 81-422-11 |
| Spare straps (pair) | 81-424-16 |



- | | |
|--|-------------|
| Lightweight leg section Plus with 50mm
mattress - For use with T Series tables only | TA-030-1077 |
| Lightweight leg section Plus with 80mm
PRS mattress - For use with T Series tables only | TA-030-1079 |
| 50mm mattress for lightweight leg
section Plus | TA-050-2200 |
| 80mm PRS mattress for lightweight leg
section Plus | TA-050-2205 |



- Douche tray complete with strainer
Compatible with T & RX series tables 81-367-85
- Douche tray for MR, J3 & J4 tables 81-367-50
- Spare strainer (fits both variants) 81-367-77



- Instrument tray for T Series tables
(Detachable tray) 81-345-06
- Instrument tray for RX500, MR and J3
(not shown) 81-342-00



- Instrument table with detachable tray
(420mm x 300mm) - Suitable for all tables 81-344-13
- Requires 1 drop handle clamp (TA-020-1084/5)



- Turner Warwick kidney supports (Pair) 81-490-03
- Requires 2 drop handle clamps (TA-020-1084/5)



- Side support for T50 (short trunk) 81-368-57
- Side support for T50 (backrest) 81-368-66



- Kidney elevator
(width adjustable, 490mm-550mm) R10-480

Also available (not shown)

- Pad for T Series lightweight leg section TA-050-1100
- Lightweight leg section for RX tables with pad 81-610-70
- Spare pad (50mm) 81-610-86
- Lightweight leg section for MR tables (bare) 81-610-78
- Pad for MR lightweight leg (50mm) 81-610-38
- Velcro straps for Direct placement legholders
(set of 3) 81-270-55




Orthopaedic and Trauma

Modular traction and positioning accessories

A range of orthopaedic accessories is available for the Eschmann T30 operating table, enabling surgeons to undertake a number of specialist procedures including leg and hand traction, tibia nailing and lateral femoral nailing.

T30 operating tables can be easily converted to enable the attachment of orthopaedic accessories using the T30 orthopaedic trunk (see page 18). A docking cart is also available to ensure quick and easy change over from a standard lower trunk to a T30 orthopaedic trunk.



T30 Orthopaedic accessory trolley

Orthopaedic Trolley - bare (UK)	TA-040-1135
Orthopaedic Trolley - bare (Euro)	TA-040-1142
Orthopaedic Trolley - bare (US)	TA-040-1181
Orthopaedic Trolley - bare (Denyer)	TA-040-1169
Orthopaedic Trolley complete with T30 orthopaedic accessories (UK)	TA-040-1136
Trolley complete with T30 orthopaedic accessories (Euro)	TA-040-1153

Complete trolley includes:

- T30 orthopaedic accessory trolley bare
- T30 pre-operative legs complete (pair)
- T30 hip sections complete (pair)
- T30 traction beam x 2
- Intermediate knuckle x 2
- Perineal post, short x 1
- Perineal post mattress short x 1
- Orthopaedic boot (pair)
- Offset traction carrier x 2
- T30 traction units x 2
- Secondary clamp x 2
- Insulated traction bracket x 1
- Tibia nailing set complete x 1
- T30 accessory slider x 2



New Product

Traction boot - quick release

Full Kit 81-943-61

Traction boot - quick release

81-943-26

Adaptor - traction boot

81-943-29

Kit Contents



Traction boot - quick release
81-943-26



Disposable traction boot liner
(pack of 6) 81-943-32



Adaptor - traction boot
81-943-29

T30 Tibia nailing kit

Full Kit TA-120-2186

Tibia nailing attachment

TA-120-2145

Tibia nailing beam

TA-110-1168

Femoral condyle pads

81-933-39

Kit Contents



TA-120-2145



81-933-39

TA-110-1168



T30 Single Leg Traction accessory kit

Full Kit TA-120-5178 - UK
TA-120-5192 - Euro
TA-120-5180 - US
TA-120-5175 - Denyer

Traction unit

TA-110-3138

Secondary clamp

TA-130-1143

Offset traction carrier

TA-110-2147

Accessory slider

TA-1104-144

Orthopaedic boot (single)

81-943-23


Intermediate knuckle

TA-110-5148

Traction beam

TA-110-1137

Kit Contents



TA-110-3138

TA-130-1143

TA-110-2147

TA-1104-144

81-943-23 (single)

TA-110-5148

TA-110-1137





Eschmann **operating**table accessories



T30 Orthopaedic system - kits and spares

T30 Infill section - perineal cut-out with 50mm pad - (UK)	TA-030-5154	Hand traction device	81-601-68
Spare 50mm pad with cut-out for T30 Infill section	TA-050-1177	Perineal post with pad (short)	TA-120-4153
T30 Infill section - standard with pad - (UK)	TA-030-5155	Spare pad for perineal post (short)	TA-050-1132
T30 Infill section bare (UK)	TA-030-5161	Perineal post with pad (long)	TA-120-4154
TA-030-5172 (Denyer) TA-030-5183 (US)		Spare pad for perineal post (long)	TA-050-1154
TA-030-5194 (Euro)		T30 Orthopaedic trunk hip sections (UK)	TA-120-1155
Spare 50mm pad for T30 Infill section	TA-050-1129	TA-120-1170 (Denyer) TA-120-1185 (US)	
Pre-operative leg supports (pair)	TA-030-4142	TA-120-1196 (Euro)	
Telescopic support	81-791-15	Accessory slider (UK)	TA-1104-144
Requires drop handle clamp (TA-020-1084/5)		TA-110-4171 (Denyer) TA-110-4182 (US)	
Inline Insulated traction bracket	TA-130-3139	TA-110-4193 (Euro)	
		Lateral femoral nailing set	TA-120-3158
		Orthopaedic boot (pair)	81-943-09

T30 Orthopaedic trunk and docking cart



T30 docking cart and orthopaedic trunk with 50mm mattress - UK/Denyer	TA-040-3134
T30 orthopaedic trunk with 50mm mattress - UK	TA-100-2001
50mm mattress for T30 orthopaedic trunk	TA-050-1128
T30 orthopaedic trunk bare	TA-100-2173 (Euro) TA-100-2181 (US)
	TA-100-2169 (Denyer) TA-100-2000 (UK)
T30 docking cart (bare) - UK/Denyer	TA-040-3133
T30 docking cart (bare) - US/Euro	TA-040-3192

T30 Carbon Traction Beams



T30 Carbon Fibre Traction Beam. UK side rails - TA-110-6159

T30 Offset Traction Carrier for Carbon Fibre Traction Beam - TA-110-6163

T30 Accessory Slider for Carbon Fibre Traction Beam. UK side rails - (Not shown) - TA-110-6174

T30 Carbon Fibre Imaging Extension



280kg

Maximum
Patient weight

Carbon Fibre Extension with 50mm mattress
(1160mm x 564mm) - TA-030-1255

Carbon Fibre Extension with 80mm PRS mattress
(1160mm x 564mm) - TA-030-1257

Replacement 50mm mattress - TA-050-6100

Replacement 80mm PRS mattress - TA-050-6103



Orthotec 5 system accessories

The Orthotec 5 system is a range of orthopaedic and traction accessories specifically developed for use with Eschmann MR, J4 and RX500 tables.

Orthotec 5 is only available for export customers.

Orthotec 5 attachment



Orthotec 5 (for use with MR & J4 tables)

81-905-18

Complete with sacral support assembly, 1 x Perineal post (81-947-26), 1 x Sacral support large (81-938-19), 1 x Sacral support extra large (81-922-19), 1 x 80mm pad (81-922-94), 1 x set of 2 Bootees, adult (81-943-09)

Orthotec 5 attachment for RX500, as 81-905-18 but with Orthotec shoe

81-908-20

Orthotec shoe for RX500

81-928-20

Orthotec 5 spares

Offset insulated traction bracket 81-930-45

Kirschner stirrup adaptor for use with 81-930-45 81-930-53

Tibia nailing accessory 81-933-04

Hand traction device 81-601-68

Sacral support plate - large 81-938-19

Sacral support plate - extra large 81-922-19

Other spares

Spare pad for counter traction bar of tibia nailing set (81-933-04) 81-933-20

Spare pad for knee stabiliser bar of tibia nailing set (81-933-04) 81-933-39




Shoulder positioners

Power assisted movement

The T20 and T30 shoulder positioners are easily attached to the respective operating tables and are suitable for all patient sizes up to 220kg. The positioners are designed to facilitate all procedures involving the shoulder girdle and upper arm. The T20 and T30 shoulder positioners utilise the powered functions of the operating table to position

the patient, minimising manual handling. Individual removable back sections ensure the patient is adequately supported while achieving unrivalled access for surgery. Fully adjustable, the head support ensures anatomically correct patient positioning.



T20 shoulder positioner

- T20 shoulder positioner 50mm mattress- 81-267-71
To fit UK/Euro side rails (clamps included)
- T20 shoulder positioner 80mm PRS mattress- 81-267-75
To fit UK/Euro side rails (clamps included)
- T20 shoulder positioner 50mm mattress- 81-268-15
To fit US side rails (clamps included)
- T20 shoulder positioner 80mm PRS mattress- 81-268-18
To fit US side rails (clamps included)
- T20 shoulder positioner 50mm mattress- 81-268-04
To fit Denyer side rails (clamps included)
- T20 shoulder positioner 80mm PRS mattress- 81-268-07
To fit Denyer side rails (clamps included)
- T20 Shoulder Positioner Upgrade Kit - 81-268-51
1 x Low Profile Head Support - Complete,
1 x 80mm Offset Head Support Column - Short (319mm)
1 x 80mm Offset Head Support Column - Long (415mm)



T30 shoulder positioner

- T30 shoulder positioner complete 50mm mattress- 81-267-64
- T30 shoulder positioner complete 80mm PRS mattress- 81-267-69
- T30 Shoulder Positioner Upgrade Kit - 81-268-54
1 x Low Profile Head Support - Complete
1 x 80mm Offset Head Support Column - Long (415mm)



Shoulder positioner spares and accessories

Replacement mattress and pad set for T20 shoulder positioner	81-533-01
Replacement mattress and pad set for T30 shoulder positioner	81-533-13
Replacement shoulder support pad for T20 and T30 shoulder positioners	81-267-73
Replacement back or head pre- support pad for T20 and T30 shoulder positioners	81-267-83
Replacement trunk and back section mattress for T20 and T30 shoulder positioners	81-267-87
Replacement head support pad for T20 and T30 shoulder positioners	81-267-96
Replacement crown or forehead strap and pad for T20 and T30 shoulder positioners	81-267-93
Standard Head Support Complete for T20 and T30 Shoulder Positioners - 1x Standard Head Support, 1x Pad for Standard Head Support, 1x Crown or Forehead Strap & Pad and 5x Chin Strap	81-268-73
Pad for Low Profile Head Support (T20 and T30 Shoulder Positioners)	81-268-57
Forehead Strap 100cm & Pad (use with low profile head support)	81-268-64
Single use strap care kit consisting - 1x Torso strap, 1x arm strap, 1x headrest support strap	81-267-32
Ultra patient care kit consisting - 1x multiflex pad set, 1x forehead strap, 1x chin strap, 1x arm strap and 1x torso strap (3 per case)	81-267-45
Storage cart for Ultra shoulder positioner	81-267-36
Securing strap for Ultra shoulder positioner	759839
Chin strap suitable for all Eschmann shoulder positioners	759753
Counter-traction strap	TA-0000-104

Low profile head support



Low Profile Head Support Complete for T20 and T30 Shoulder Positioners - **81-268-61**
Kit contains:
 1 x Low Profile Head Support, 1 x Pad for Low Profile Head Support, 1 x Forehead Strap & Pad and 5 x Chin Strap - (not shown)

Ultra shoulder positioner

Shoulder positioning support

81-267-13

Requires blade clamp:

UK/Euro sidebar **81-273-49**

US Sidebar **81-273-52**



180kg

Maximum
Patient weight



Orthopaedic limb accessories



- Weightless Shoulder Suspension System including rotary clamp (UK) 81-274-02
- Weightless Shoulder Suspension System including rotary clamp (Euro) 81-274-05
- Weightless Shoulder Suspension System including rotary clamp (US) 81-274-07
- Sterile Arm Traps (case of 6) for use with Weightless Shoulder Suspension System (81-274-02) 81-274-13
- Sterile Impervious Arm Traps (case of 6) for use with Weightless Shoulder Suspension System (81-274-02) 81-274-24




Counter traction brace for beachchair/shoulder surgery positioning
Requires 1 rotary anti drift clamp (TA-020-1082/6/7) TA-090-1109



Prepper post, for temporary limb positioning (775mm post)
Requires 1 rotary accessory clamp (81-270-45/48/52) 81-270-30



Arthroscopic Stress Post 81-271-30
Requires 1 rotary accessory clamp (81-270-45/48/52)



Arthroscopic legholder 81-270-95
Accepts legs up to 890mm in circumference
Requires 1 rotary accessory clamp (81-270-45/48/52)

Spinal and Neurosurgery

Versatile, safe and secure


The range of Doro® neurosurgical accessories from Eschmann allows safe, secure patient positioning for a wide variety of neurosurgery procedures.



Standard system - base units

DORO® Adjustable base unit Requires relevant accessory attachment bar to fit table (see below)	81-590-01
DORO® Adjustable base unit for Eschmann tables (Not shown) Fits directly onto RX, J3 and MR models	81-590-12
DORO® Adjustable base unit for T Series tables Fits directly onto T Series tables	81-594-04
Accessory Attachment bar For use with T20/30 to allow narrow pin attachments	TA-030-3004
Accessory Attachment bar For RX system, MR and J3	81-574-05


81-590-01



81-594-04



TA-030-3004





Standard system

DORO® Crossbar adaptor For procedures requiring seated positions, adjustable. Requires Doro® adjustable base unit (81-590-01) and 2 rotary accessory clamps 81-270-45 (UK) 81-270-52 (Euro) 81-270-48 (US)	81-593-20
DORO® 6" Transitional Member	81-590-23
DORO® Skull Clamp	81-590-56
DORO® Swivel adaptor	81-590-34
DORO® Ball and Socket swivel adaptor	81-590-45
DORO® Horseshoe Headrest (Adult) with extension bar and gel pads	81-591-11
DORO® Horseshoe headrest (Paediatric) with extension bar and gel pads	81-591-44
DORO® Horseshoe pad, Adult, left	81-591-22
DORO® Horseshoe pad, Adult, right	81-591-33
DORO® Horseshoe pad, Paediatric, left	81-591-55
DORO® Horseshoe pad, Paediatric, right	81-591-66
DORO® Gel Prone headrest	81-591-77



81-590-45



81-590-56



81-590-67




81-590-69

Multi purpose system

DORO® Multi purpose Skull Clamp, Adult and Paediatric (set) comprising: Multi purpose skull clamp, Head support, Pin holder paediatric, Single pin holder, Pin holder adjustable - paediatric, Disposable skull pins - paediatric (pack of 36), Trio gel pad adjustable, Gel head rings with base 90mm (3 pcs), Gel ear ring with base 90mm (2 pcs), Gel head ring with base 140mm (1 pc), Adjustable key.	81-590-67
DORO® Gel head ring with base, 90mm	81-590-69
DORO® Gel head ring with base, 140mm	81-590-71

Radiolucent system

DORO® Radiolucent headrest complete system Constructed of Novotex®. No metal parts to interfere with intra-operative imaging. Set comprising of: Storage case, base unit, transitional member-long, transitional member-short Swivel adaptor, skull clamp, disposable titanium pins - Adult (36 pieces).	81-591-88
DORO® Transitional member Aluminium base unit	81-592-21
DORO® Transitional member Radiolucent, long (160mm)	81-592-32
DORO® Transitional member Radiolucent, short (120mm)	81-592-43
DORO® Radiolucent swivel adaptor	81-592-54
DORO® Radiolucent interface for skull clamp	81-592-65
DORO® Aluminium interface for skull clamp	81-592-76
DORO® Swivel Horseshoe headrest, Radiolucent, Adult, with extension bar	81-591-80
DORO® Swivel Horseshoe headrest, Radiolucent, Adult, without extension bar	81-591-82
DORO® Swivel Horseshoe headrest, Radiolucent, Paediatric, with extension bar	81-591-84
DORO® Swivel Horseshoe headrest, Radiolucent, Paediatric, without extension bar	81-591-86
DORO® Radiolucent skull clamp	81-592-87




Autoclavable System

DORO® Teflon adjustable base unit Made from aluminium with Teflon coating. Requires accessory attachment bar (TA-030-3004) For use with T Series tables.	81-593-42
DORO® Teflon swivel adaptor	81-593-53
DORO® Teflon skull clamp	81-593-64



Skull pin sets

DORO® Adult reusable skull pins (set of 3)	81-590-78
DORO® Paediatric reusable skull pins (set of 3)	81-590-89
DORO® Disposable adult skull pins Colour coded blue for ease of use, sterile packs of 3 pins (36 pins/box)	81-591-00
DORO® Disposable paediatric skull pins colour coded yellow for ease of use, sterile packs of 3 pins (36 pins/box)	81-592-98
DORO® Disposable adult titanium skull pins colour coded turquoise for ease of use, sterile packs of 3 pins (36 pins/box)	81-593-09





Spinal Frames



- | | |
|---|-----------|
| Radiolucent Spine Surgery Frame System
(includes Cart 81-271-48, Crank Arm 81-271-45
and case of 6 sets Skin Care Covers 81-271-36) | 81-271-33 |
| Skin Care Covers (case of 6 sets) | 81-271-36 |
| Standard Covers (case of 6 sets) | 81-271-39 |
| Replacement Pad Set | 81-271-42 |
| Crank Arm for Radiolucent Spine Surgery Frame | 81-271-45 |



- | | |
|--|-----------|
| Laminectomy support pad (montreal mattress)
(990mm x 130mm) | 81-569-13 |
|--|-----------|



- | | |
|--|-----------|
| Wilson spinal surgery frame
Complete with ProneView® head support system | 81-569-58 |
| Patient care kit - single patient use (6 per case)
Kit comprises:
1 x ProneView® cushion insert, 2 x pad covers, 1 x
base cover drape, 2 x foam arm cradles, 1 x cross
bar adaptor | 81-569-66 |
| ProneView® adjustable mirror | 81-569-04 |
| Replacement ProneView® Helmet | 81-569-47 |
| Replacement Gel pad set | 81-569-25 |
| Replacement Crank handle with torque
limiter | 81-569-75 |

Patient Positioning Aids

Flexible patient positioning


The Modular Patient Positioning System (MPPS) from Eschmann has been designed to complement the range of surgical procedures carried out in today's operating theatres by utilising a range of modular components that can be easily configured to hold and position the patient in tune with their anatomy.



Simple, fast, configuration

The MPPS system has been designed with speed and ease of use in mind. Simply changing the support pad and a small adjustment of the support bar can provide outstanding support for a completely different position or procedure.

Make up your own MPPS set or choose one of the MPPS kits overleaf.



Support bar (280mm)	TA-020-1040
Support bar (440mm)	TA-020-1069
Also available, not shown	
Support bar left hand side (600mm)	TA-020-1083
Support bar right hand side (600mm)	TA-020-1068



Ball locking clamp assembly (150mm)	TA-020-1039
Ball locking clamp assembly (250mm)	TA-020-1070
Ball locking clamp only	TA-070-2071



Small circular support
(100mm Diameter)

TA-050-1041

Large circular support
(150mm Diameter)

TA-050-1088



Small rectangular (220mm)

TA-050-1042

Large rectangular (300mm)

TA-050-1072



Square support (180mm x 210mm)

TA-050-1061



Knee support (single) (bare)

TA-070-1075

Spare pad for Knee support (Standard)

TA-050-1066

Spare pad for Knee support (large)

TA-050-1067



Small armboard (400mm)

TA-050-1059

Large armboard (575mm)

TA-050-1058



Carter Braine style forearm support (bare)

TA-070-1062

Spare pad for Carter Braine style
forearm support

TA-050-1063



Upper arm support (bare)

TA-070-1064

Spare pad for upper arm support

TA-050-1065

MPPS Kits

Arm support set




Arm support set (575mm)
TA-020-1048

Consisting of 575mm armboard (TA-050-1058)
Support bar assembly (TA-020-1060)
Retaining strap (TA-000-0073)
Set requires 1 rotary anti drift clamp (TA-020-1082/6/7)

Arm support set (400mm)
TA-020-1049

As above with 400mm armboard

Chest/Pelvis square support set



Chest/pelvis square support set
TA-020-1050

Consisting of square chest/pelvis support pad
180mm x 210mm (TA-050-1061)
Support bar 280mm (TA-020-1040)
Ball locking clamp assembly 150mm (TA-020-1039)
Set requires 1 drop handle clamp (TA-020-1084/5)

Carter Braine style forearm support set



Carter Braine style
forearm support set
TA-020-1051

Consisting of forearm support (TA-070-1062)
Forearm support pad (TA-050-1063)
Support bar 280mm (TA-020-1040)
Retaining strap (TA-000-0073)
Set requires either 1 drop handle clamp (TA-020-1084/5) or
1 rotary anti drift clamp (TA-020-1082/6/7)

Upper arm support set



Upper arm support set
TA-020-1052

Consisting of upper arm support (TA-070-1064)
Upper arm support pad (TA-050-1065)
Support bar 280mm (TA-020-1040)
Retaining strap (TA-000-0073)
Set requires either 1 drop handle clamp (TA-020-1084/5) or
1 rotary anti drift clamp (TA-020-1082/6/7)

Adjustable pelvis support set



Adjustable pelvis support set
260mm support arms TA-020-1053

Consisting of 2 small circular support pads (TA-050-1041)
Pelvis support bar assembly (TA-020-1074)
Set requires 1 drop handle clamp (TA-020-1084/5)

Adjustable pelvis support set
390mm support arms TA-020-1026

Extra long support set as above with 390mm support arms



Knee support set



Standard knee support set (Pair)
TA-020-1054


Large knee support set (Pair)
TA-020-1055

Consisting of 2 knee supports (TA-070-1075)
2 x standard knee support pads (TA-050-1066)
Support bar left hand side 600mm (TA-020-1068)
Support bar right hand side 600mm (TA-020-1083)
2 x retaining straps (TA-000-0073)

As above with 2 x large knee support pads (TA-050-1067)

Sets require 2 rotary anti drift clamps (TA-020-1082/6/7)

Rectangular lateral support set



Rectangular lateral support set
TA-020-1056

Consisting of rectangular support 220mm (TA-050-1042)
Support bar 280mm (TA-020-1040)
Ball locking clamp assembly 150mm (TA-020-1039)
Set requires 1 drop handle clamp (TA-020-1084/5)

Circular lateral support set



Circular lateral support set
TA-020-1057

Consisting of circular support pad (TA-050-1041)
Support bar 280mm (TA-020-1040)
Ball locking clamp assembly 150mm (TA-020-1039)
Set requires 1 drop handle clamp (TA-020-1084/5)

Circular chest/pelvis support set



Circular chest/pelvis support set
TA-020-1091

Consisting of large circular support pad 150mm
(TA-050-1088)
Support bar 280mm (TA-020-1040)
Ball locking clamp assembly 150mm (TA-020-1039)
Set requires 1 drop handle clamp (TA-020-1084/5)

Gel pads

Patient positioning comfort

Eschmann offers a wide range of pressure relieving gel pads for use in patient support and positioning. All gel pads are made from 100% silicone based gel, which resembles the texture of human tissue, providing exceptional comfort and superior performance.



Contoured armboard gel pad
(460mm x 180mm x 35mm) **81-302-72**



Perspex arm support gel pad
(170mm x 110mm x 15mm) **81-368-31**




Patient positioning gel strap
(Max length approx. 2500mm) **81-462-24**



Sacral support gel pad
(450mm x 445mm x 30mm) **81-469-17**



Gel Head stabilisers
Adult (260mm x 250mm x 50mm) **81-374-59**
Child (150mm x 35mm) **81-374-78**
Paediatric (90mm x 20mm) **81-374-89**



Universal gel pads
Small (160mm x 80mm x 45mm) **81-379-87**
Medium/small (330mm x 70mm x 50mm) **81-379-75**
Medium/large (300mm x 120mm x 60mm) **81-379-63**
Large (360mm x 125mm x 70mm) **81-470-39**



Eschmann **operating**table accessories



'D' shaped gel positioner pad, general purpose (350mm x 150mm x 75mm)

81-469-15



Prone headrest gel pad (280mm x 240mm x 140mm)

81-469-13



Operating table gel pad - full length (1800mm x 520mm x 10mm)

81-469-39

Operating table gel pad - short (520mm x 520mm x 10mm)

81-469-63

Operating table gel pad 3/4 length (1150mm x 520mm x 10mm)

81-469-75



Operating table gel pad, 3/4 length with perineal cutout (1150mm x 520mm x 10mm)

81-469-54



Contoured leg protector (515mm x 180mm x 40mm)

81-470-61



Retaining strap (Max length 800mm)

TA-000-0073



Heel protector (190mm x 175mm x 30mm)

81-470-45



Perineal post gel pad (300mm x 260mm x 10mm)

81-469-42

Replacement handsets

Replacement handsets and mains cables

The following range of replacement corded handsets are available. Please ensure that the correct handset is selected for the correct type of operating table and sidebar.

For replacement handsets for RX range please contact Eschmann.

Mains cables

4 metre 3 core mains cable with 3 pin UK plug	714188
4 metre mains cable with Euro plug (not fused)	715254
4 metre mains cable with Australian mains plug (not fused)	391177
2 metre mains cable with North American 3 pin plug	111589

Table handsets

T50 Models	
Handset with UK sidebar hook	T00-570-0012
Handset with Euro/US sidebar hook	T00-570-0024
T30 Models (all variants)	
Handset with UK sidebar hook	T00-563-0068
Handset with Euro/US sidebar hook	T00-563-0079
T20 Models (all variants)	
Handset with UK sidebar hook	T00-562-0024
Handset with Euro/US sidebar hook	T00-562-0035
T10 Models (all variants)	
Handset with UK sidebar hook	T00-564-0103
Handset with Euro/US sidebar hook	T00-564-0114





T Series operating tables

Eschmann operating tables have been trusted by medical professionals the world over for over 100 years. The T Series family of operating tables is no exception, providing outstanding choice and flexibility for both general and specialist surgical procedures.

Basic General Surgery Operating Table

- Capable of 250kg maximum patient weight
- 703-1003mm height range (excluding mattress)
- 600mm wide table top
- Full length radiolucent top with integral X-ray cassette tunnel
- Interchangeable head and leg sections
- Lightweight, seamless mattresses with secure "Velcro free" fixings
- Choice of sidebars available - maximising use of existing accessories



T10



General Surgery Operating Table

- Capable of 450kg maximum patient weight
- 670-1070mm height range (excluding mattress)
- 600mm wide top plus optional width extenders
- Full length radiolucent top with integral X-ray cassette tunnel
- Imaging made easy with 250mm traversing top
- Interchangeable head and leg sections
- Robust corded handset with battery status indicator
- Lightweight, seamless mattresses with secure "Velcro free" fixings




T20



Find further information on T
Series operating tables
www.eschmann.co.uk

General Purpose & Orthopaedic Modular Table System

- Capable of 450kg maximum patient weight
- 670-1070mm height range (excluding mattress)
- 600mm wide top plus optional width extenders
- Simple to configure for general orthopaedic surgery
- Wide range of orthopaedic attachments available
- Compact accessory trolley and docking cart for trunk sections - reduces manual handling
- Soft start/stop technology
- Choice of sidebars available - maximising use of existing accessories
- Shoulder surgery made easy with T30 shoulder positioner



Obstetrics and Gynaecology Operating Table

- Shorter table top means less manual handling
- 300kg maximum patient weight in lithotomy position
- Powered back rest allows for 'Assisted Delivery' position
- Increased speed on auto level function
- Weight assisted articulating leg sections
- Improved surgical access via base cut-out





General Information

Ordering

Please ensure you order the correct clamps for accessories that require attachment to the table sidebars. If you need help finding out if your existing clamps will be suitable please contact your local Eschmann representative who will be happy to help.

Service

Eschmann products are supported by a worldwide network of fully trained engineers, offering high quality Eschmann spare parts. For further information on the range of Service Contracts available, please contact your local Eschmann representative.

In the UK all Eschmann operating tables are fully supported by our team of professional engineers.



Find further information on
new product releases
www.eschmann.co.uk

Imagery in this publication has been taken to illustrate the operating table movements and as such do not constitute a guide to correct patient positioning

All accessories supplied are CE marked in compliance with the Medical Devices Directive 93/42/EEC

Products are not made with natural rubber latex

Distributor Information

Global Leaders in Healthcare



ESCHMANN