



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Baxter Sigma Spectrum Infusion Pump - Refurbished



- Intuitive Design
- Bright LCD Display
- Small/Compact Design
- Standard IV Sets
- Multiple Safety Systems
- Audible Alarm
- Built-In Drug Library
- Rechargeable Battery
- Completely Refurbished

The Baxter Sigma Spectrum Infusion Pump provides rapid, intuitive infusions in a small package. The unit uses standard administrations sets for both gravity and pump infusions, lowering costs. Smart technology in order to help prevent pump related medical errors. The bright and easy to read screen is easy to read and also display on-board help screens. The Baxter Sigma Spectrum also has wireless options and is EMR compatible.

The Baxter Sigma Spectrum comes completely refurbished with battery, pole clamp, power cord, users manual, and a comprehensive warranty.

Specifications:

- Dimensions: 64.2"D x 4.2"W x 6.3"H
- Volume Range: 0.1 to 99.9 mL (in 0.1 mL increments)
- Volume Range: 100 to 999 ML (in 1.0 mL increments)
- Battery: LiON
- Weight: Approx 2 lbs
- Mounting Type: Pole Clamp
- Power Requirements:
 - 120 VAC, 60Hz

VENTUREMEDICAL.COM

© 2007-2018 Venture Medical ReQuip, Inc.