



www.venturemedical.com

Venture Medical ReQuip, Inc.  
6008 Bonacker Drive  
Tampa, FL 33610  
800-627-3215 | 813-740-8640  
Fax: 800-466-1251  
sales@venturemedical.com

## Cardiac Science Powerheart G3 AED - Refurbished



- Biphasic
- Intuitive Voice Prompts
- Lightweight
- Rescue Ready Visual Indicator
- A One-piece Electrode Pad for Fast and Accurate Placement
- Lithium-Ion Battery for Multiple Events
- Extremely Portable
- Built in Carry Handle
- Pad Storage Inside Case

The Powerheart AED is a biphasic, user-friendly, advanced AED that features one button operation and Rescue Ready technology that allows clinicians to know that the AED is ready for use.

The Rescue Coach voice prompts the user through the process in an emergency. The bold, authoritative voice is aimed at keeping the caregiver calm and focused, while being audible enough to hear in a potentially loud rescue environment. The voice prompts means more confidence and less confusion for the caregiver.

The Cardiac Science Powerheart G3 comes complete with new battery, set of adult pads, users manual.

The Powerheart G3 comes completely refurbished with a comprehensive warranty.

### Specifications:

- Dimensions: 3.3"h x 10.6"w x 12.4"d
- Energy Range: 150J -360J, Escalating
- CPR: Metronome for Compression Frequency
- Charge Hold Time: 30 seconds
- Waveform: Rectilinear Biphasic
- Weight: 6.6 lbs
- Power: Lithium Ion Battery