

## Valleylab LigaSure Vessel Sealing Unit Refurbished



- Digital Display
- Bipolar Mode
- Macrobipolar Mode
- Constant, Consistent Output
- Average Seal Cycle is 2-4 Seconds
- Foot-Switch Included
- Completely Refurbished

The LigaSure™ vessel sealing system gives you something no other electro-surgical tool can - permanent vessel occlusion. The LigaSure™ system's revolutionary technology replaces almost all other hemostatic tools because it doesn't just ligate vessels, it actually fuses vessel walls to create a permanent seal.

Included with the unit: Footswitch, Operators Manual and Power Cord. The Valleylab Ligasure comes fully refurbished with a comprehensive warranty.

### Specifications:

- Macrobipolar: 473 kHz sinusoid
- Bipolar: 473 kHz Sinusoid
- Seal Mode: 473 kHz Sinusoid, Pulsed
- Height (unit only): 12.7 cm (5")
- Width (unit only): 38.6 cm (15.2")
- Length (unit only): 40.6 cm (16")
- Weight (unit only): 5.9 kg (13 lbs)
- Power Requirements: 120V