

UltraScan Versa™ Premier

The ultimate multi-modality ultrasound platform

MEDICAL POSITIONING, INC.
www.medicalpositioning.com
800-593-3246



STANDARD FEATURES

| | | | | | |
|------------------------------|---|--------------------------|---|---|----------------------|
| Height Adjustability: | Maximum Patient Capacity: | Fowler (Head Up): | Trendelenburg/Reverse Trendelenburg: | Calf Section: | Braking: |
| 22" to 38" | 500 lbs | 0 - 90° | 15° /25° | 0 - 83° | Single Pedal Braking |
| Control System: | Foot Down Section: | Stirrups: | Handrails: | Arm Boards: | |
| Hand Controller | Electrically adjustable Built-In Folding Foot Platform | Adjustable and Hideaway | Removable with Hideaway Bracket | Dual Articulating Patient Arm Boards and Back Support | |

PRODUCT SPECIFICATIONS

| | | |
|-----------------------|-----------------------------------|---|
| Length: | Width: | Weight: |
| 72.5" | 28.5" | 350 lbs |
| Foam: | Electrical: | Vinyl: |
| TB 117-2013 Section E | 100 - 240 VAC, 50/60 Hz, 4.5A Max | TB 117-2013 NFPA 260 Cover Fabric - Class 1 UFAC Fabric - Class 1 |

PRODUCT BENEFITS

Perform multiple ultrasound procedures using one table


Optimize patient access and positioning to help reduce image acquisition time

Help protect sonographers from potential career-ending injuries by improving ergonomics

Bariatric capable

Increase lower extremity vessel dilation

CERTIFICATIONS

| | | |
|--|------------------|---|
| Industry Standards: | Warranty: | Standard Color: |
| UL ANSI/AAMI ES60601-1 (2005) + AMD (2012), CAN/CSA C22.2 No. 60601 (2014) | 5 Years |  |

OPTIONAL FEATURES

Carotid Head Rest

Foot Control

Head Rail

IV Pole and Holder

Paper Roll Holder and Cutter

Push Handles

Safety Straps

Underbed LED Lighting