


# Stress EchoBed<sup>®</sup>

The ultimate platform for Supine Bike Stress Echocardiography

MEDICAL POSITIONING, INC.  
www.medicalpositioning.com  
800-593-3246




## STANDARD FEATURES

| Height Adjustability:                   | Maximum Patient Capacity:  | Drop Sections:                                      | Trendelenburg/Reverse Trendelenburg:   | Lateral Tilt Left:  |
|---|----------------------------|---|--|---|
| 24" to 34"                              | 500 lbs                    | Imaging Drop Section<br>Sonographer Drop<br>Section | 15°/15°  | 40°   |
| Control System:                         | Braking:                   | Ergonomic Features:                                 | Additional Features:   |   |
| Hand Controller with Memory Positioning | Individual Locking Casters | 2-Way Sonographer Drop Section/Patient Back Support | Safety Rails<br>Safe-T-Wedge<br>Patient Safety Strap<br>Paper Roll Holder and Cutter<br>Shoulder support | Removable Ergometer with computer controller designed for stress echocardiograms, including 10 preset protocols, 4 programmable protocols, 1 manual protocol, and a patient coach with visual and auditory prompts. |

## PRODUCT SPECIFICATIONS

| Length:                                       | Width:                      | Weight:   |
|---|-----------------------------|---|
| 78"   | 30"                         | 465 lbs   |
| Foam:   | Electrical:                 | Vinyl:  |
| California Technical Bulletin 117 - Section E | 120 VAC, 50/60 Hz, 2.6A Max | California T.B. 117-2013<br>NFPA 260- Cover Fabric - Class I<br>UFAC Fabric - Class I |

## CERTIFICATIONS

| Industry Standards:                 | Warranty: | Standard Color:   |
|-------------------------------------|-----------|---|
| UL 60601-1, CAN/CSA C22.2 No. 601.1 | 5 Years   |  |

## PRODUCT BENEFITS

|   |  |
|---|--|
| Induces symptoms at lower heart rates   | Avoids the need for post-exercise transfer                       |
| Provides imaging during intermediate stages to identify the ischemic threshold          | Evaluates diastolic pressures and manages heart failure patients |
| Improves diagnostic capabilities for valvular pathology and pulmonary vascular problems | 4 point safety system  |
| Eliminates rapidly falling heart rates seen post-treadmill                              | Bariatric capable  |

## OPTIONAL FEATURES

|  |                       |
|--|-----------------------|
| Ergometer Storage Cart                           | IV Pole and Holder    |
| Ergometer Insert Pediatric/<br>Geriatric Adapter | Padded Arm Board      |
|  | Sonographer Extension |