

UltraScan Versa™ Table

The ultimate ultrasound platform

MEDICAL POSITIONING, INC.
www.medicalpositioning.com
800-593-3246

Weight capacity up to 500 lbs

90° Fowler positioning

Optional +15°/-15° Trendelenburg
or +15°/-25° Trendelenburg
positioning

Adjustable hideaway stirrups

22" to 38" height

Electric built-in calf section
and optional folding foot
platform



ergonomic * innovative * proven

Optional features may be shown. Actual color may vary depending on monitor resolution and print quality.

STANDARD FEATURES

Height Adjustability:	Maximum Patient Capacity:	Fowler (Head Up):	Calf Section:	Braking:
22" to 38"	500 lbs	0 - 90°	0 - 83°	Individual Locking Casters
Control System:	Foot Down Section:	Stirrups:		
Hand Controller	Electrically Adjustable Optional Built-In Folding Foot Platform	Adjustable and Hideaway		

PRODUCT SPECIFICATIONS

Length:	Width:	Weight:
72.5"	28.5"	350 lbs
Foam:	Electrical:	Vinyl:
TB 117-2013 Section E	100 - 240 VAC, 50/60 Hz, 4.5A Max	TB 117-2013 NFPA 260 Cover Fabric - Class 1 UFAC Fabric - Class 1

PRODUCT BENEFITS


Perform multiple ultrasound procedures using one table

Optimize patient access and positioning to help reduce image acquisition time

Help protect sonographers from potential career-ending injuries by improving ergonomics

Bariatric capable

CERTIFICATIONS

Industry Standards:	Warranty:	Standard Color:
UL ANSI/AAMI ES60601-1 (2005) + AMD (2012), CAN/CSA C22.2 No. 60601 (2014)	2 Years	

OPTIONAL FEATURES

+15°/-15° Trendelenburg Positioning

+15°/-25° Trendelenburg Positioning

Handrails with Hideaway Brackets

IV Pole and Holder

Paper Roll Holder and Cutter

Underbed LED Lighting Foot Control

Foot Board

Single Pedal Braking

Push Handles

Carotid Head Rest

Head Rail

Safety Straps

Battery Functionality