




The ultimate platform for video swallow studies

MEDICAL POSITIONING, INC.  
www.medicalpositioning.com  
800-593-3246



Optional features may be shown. Actual color may vary depending on monitor resolution and print quality.

## STANDARD FEATURES

<b>Height Adjustability:</b>	<b>Maximum Patient Capacity:</b>	<b>Fowler:</b>	<b>Calf Section:</b>	<b>Arm Rests:</b>	<b>Radiolucent Imaging:</b>	<b>Seat Rotation:</b>
19" to 33"	600 lbs	90°	0° to 75°	Fold Away and Height Adjustable	32" x 13"	90° Left and Right
<b>Powered Operation:</b>	<b>Control System:</b>	<b>Braking:</b>	<b>Additional Features:</b>			
Battery Operated	Hand Controller	Single Pedal Braking	Oxygen Tank Holder Patient Safety Straps Transport Handles Fold away seat back Fold away and height adjustable arm rests			

## PRODUCT SPECIFICATIONS

<b>Table Length:</b>	<b>Chair Length:</b>	<b>Width:</b>	<b>Weight:</b>
79.25"	33"	26"	365 lbs
<b>Foam:</b>	<b>Electrical:</b>	<b>Vinyl:</b>	
California Technical Bulletin 117 - Section E	100-230 VAC. 50/60 Hz, 5.0A Max	California T.B. 117-2013 NFPA 260- Cover Fabric - Class I UFAC Fabric - Class I	

## OPTIONAL FEATURES

- IV Pole and Holder
- Charger and Second Battery

## CERTIFICATIONS

<b>Industry Standards:</b>	<b>Warranty:</b>	<b>Standard Color:</b>
UL ANSI/AAMI ES60601-1 (2005) + AMD (2012), CAN/CSA C22.2 No. 60601 (2014)	1 Year	 Latch GUA-1304