

■ SELECT COLLECTION

Verō


Available: Power, Manual, Power XL, Manual XL



Champion is leading the way when it comes to engineering medical seating that is both innovative and functional.

The Verō features a modular design that allows you to choose the exact configuration required to meet your specific needs.

Ergonomic design can prevent repetitive strain injuries among the medical team and caregivers while providing safety and comfort for the patient.



Key Benefits

- Multiple recline positions promote patient independence.
- The controlled recline descent provides gentle patient positioning while avoiding injury to caregivers.
- Modular design allows choices in how the chair is anchored or moved. Whether stationary or for patient transport, Verō has a perfect solution.
- For quality engineering and patient-centric innovation that prioritizes functionality, safety, and comfort, Champion products comply with ISO 13485 and are made in the USA.

Common Applications:

General Patient Room Seating, Dialysis Settings

Standard Features



Trendelenburg Positioning



Fixed or Swing Arms




5" (12.7 cm) Casters
3" (7.62 cm) Dual Casters
Faux Feet/Pedestal Feet




Foot Tray

Space Requirements


Front View




Top View



Side View



Top View with Swing Away Arms



Specifications

	Power (6P)	Power XL (6P)	Manual (6Y)	Manual XL (6Y)
Weight Capacity	350 lbs. / 158 kg	500 lbs. / 226 kg	350 lbs. / 158 kg	500 lbs. / 226 kg
Back Height (from floor ♦ from seat)	47.50" ♦ 29" 120.65 cm ♦ 73.66 cm	47.50" ♦ 29" 120.70 cm ♦ 73.70 cm	48.25" ♦ 29" 122.56 cm ♦ 73.66 cm	48.25" ♦ 29" 122.56 cm ♦ 73.70 cm
Overall Width (no tables & arms closed)	29.50" / 74.93 cm	34" / 86.36 cm	29.50" / 74.93 cm	34" / 86.36 cm
Seat Width	22.50" / 57.15 cm	26.50" / 67.31 cm	22.50" / 57.15 cm	26.50" / 67.31 cm
Seat Height	23.75" / 60.33 cm	23.75" / 60.33 cm	22.75" / 57.79 cm	22.75" / 57.79 cm
Arm Height (highest point from floor)	29" / 73.70 cm			
Seat Depth	23.50" / 59.70 cm	23.50" / 59.70 cm	22" / 55.88 cm	22" / 55.88 cm
Overall Upright Depth	41" / 104.14 cm	41" / 104.14 cm	39.50" / 100.33 cm	39.50" / 100.33 cm
Trendelenburg Depth	69" / 175.26 cm	69" / 175.30 cm	69" / 175.26 cm	69" / 175.30 cm

(Specs shown for models with 5" casters. Pedestal feet and 3" caster options available. Talk to a Sales Rep for details.)

Available Options



Black Arm Caps



Heat + Massage Interior Arm Switch



Side Tables



Pillow for Head Or Lumbar Support



IV Pole Options

(Additional available options are also offered. Talk to a Sales Rep for details.)