



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Mortara ELI 250 EKG Unit - Refurbished



- 12 Lead EKG
- Interpretation
- Full Page Printout
- LCD Display with Backlight
- Preview Screen (Waveform Display)
- Full Keyboard
- Patient Annotation
- Memory (60 Resting ECGs)
- Portable and Lightweight
- "One-Touch" ECG and Rhythm Acquisition
- EMR Compatible

The compact and lightweight Mortara ELI 250 provides a full-featured electrocardiograph.

Features include simultaneous 12 lead ECG acquisitions, "one touch" automatic ECG, interpretation, a backlit preview screen with real-time waveforms and much more. The full alphanumeric keyboard provides the ability to perform patient annotation and the internal memory can store up to 60 digital ECG records.

The Mortara ELI 250 also is EMR compatible and able to be used with Mortara's E-Scribe™ and Athena products as well as third party EMR, HIS, and PACS systems via XML, PDF, UNIPRO, and DICOM®.

Includes: Patient Cable, Leadwires, Pack of Paper, New Battery, Power Cord and Operators Manual.

The Mortara ELI 250 EKG comes completely refurbished with a comprehensive warranty.

Specifications:

- Display Type: LCD with Backlight
- Dims: 15.5" x 17" x 4"
- Weight (unit only): 11.25 lbs
- Power Requirement: 120v
- Printer Type: Thermal Array