



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Midmark 230 Power Exam Table Refurbished



- +/- Height
- Tilt
- +/- Backrest
- +/- Leg Rest
- Table: Up/Down (22.5" - 40")
- Tilt (0° - 30°)
- Foot: Up/Down (0° - 90°)
- Back: Up/Down (0° - 85°)
- Adjustable Headrest (Manual)
- Hideaway Adjustable Stirrups
- Pull-Out Cysto Pan
- Table Paper Holder
- Table Paper Strap
- Foot Switch or Hand Control
- Accessory Side Rails

The Midmark 230 delivers reliable performance and ergonomic design, making it a trusted choice for efficient outpatient procedures. Refurbished by Venture Medical, it meets high-quality standards while offering excellent value.

This power chair features height adjustment from 22.5" to 40", controlled by footswitch or hand control. It includes an articulating headrest, accessory side rail, paper holder, and a pull-out drain pan. The foot section locks into four proctology positions, with built-in heel stirrups for added versatility.

Ideal for family practice, OB/GYN, dermatology, internal medicine, and surgical specialties, the Midmark 230 enhances both patient comfort and clinical workflow.

Specifications:

- Minimum Height: 22.5" (57.2 cm)
- Maximum Height: 40" (101.6 cm)
- Table Length (foot extended): 72" (185.4 cm)
- Table Length (foot and headrest extended): 83" (210.8 cm)
- Width (seat): 28" (71.1 cm)
- Width (back): 28" (71.1 cm)
- Width (headrest): 18" (45.7 cm)
- Width (foot support): 18" (45.7 cm)
- Patient Weight Capacity: 450 lbs. (204.1 kg)
- Paper Roll Holder: Holds 18" x 3.5" paper roll (45.7 cm - 8.9 cm)
- Foot or Hand Control: Height, tilt, back, foot/leg, Stop, and auto return
- Electrical: Domestic: 115 VAC, 60Hz, 10AMP, UL Listed, CSA and CE Certified

VENTUREMEDICAL.COM

© 2007-2025 Venture Medical ReQuip, Inc.