

Midmark 625 Power Exam Table Refurbished



Stock Photo – Options May Vary

The Midmark 625 Barrier-Free power exam table is the most patient accessible exam table in the industry. The table lowers to a low height of 18" and lays flat for easy and safe patient entry and exit. The table features pull out stirrups, power height and backrest adjustment, extendable footrest with pad, hospital grade power outlet and a paper roll holder. Optional features include pelvic tilt bar and drawer warmer.

The Midmark 625 also boasts a 650 lbs patient weight capacity for bariatric patients.

The Midmark 625 comes completely refurbished with new upholstery (your choice of colors), hand control, table paper roll holder, table paper strap and a comprehensive warranty.

- Maximum Patient Accessibility
- Power height and Back
- Hand Control
- Low 18" Minimum Height
- Bariatric Weight Capacity (650 lbs)
- Pull Out Stirrups
- Extendable Leg Section with Pad
- Paper Roll Holder
- Drawer Warmer (optional)
- Pelvic Tilt (optional)
- Table Paper Holder
- Table Paper Strap
- Completely Refurbished

Specifications:

- Weight Capacity: 650 lbs
- Weight of Table: 433 lbs
- Tabletop Height: 18"- 37"
- Pelvic Tilt: 7°
- Dimensions:
 - Upholstered Top: 28" W x 59.5" L
 - With Footrest Extended: 75.5" L
- Paper Rolls: up to 21" x 3.5"
- Power: 115 VACS, 5 A Max

***Brochure and images are a general overview of the product. Actual product options may vary.*