



www.venturemedical.com

Venture Medical ReQuip, Inc.  
6008 Bonacker Drive  
Tampa, FL 33610  
800-627-3215 | 813-740-8640  
Fax: 800-466-1251  
sales@venturemedical.com

## Midmark 647 Power Podiatry Chair/Table - Refurbished



Stock Photo – Options May Vary

The Midmark 647 is built to meet the demands of today's podiatry and lower extremity care. With smooth powered adjustments, ergonomic access, and a patient-focused design, it enhances comfort while improving workflow for providers.

The Midmark 647 features a powered height system adjusts from 21" to 40", making patient transfers easier and positioning more ergonomic. A wide, flat seat surface and chair arms that rotate completely out of the way further support patient accessibility and safety.

The 647 has an expertly designed foot section with touch controls on both sides and one-handed release handles. A sliding stainless-steel debris pan adds convenience and hygiene at the point of care. The chair offers full-powered multi-positioning—adjusting height, seat, back, and tilt—supporting back angles from 0° to 80° and seat tilt up to 30°.

- Powered height adjustment raises and lowers the seat from 21" to 40"
- Expertly designed foot section with touch controls located on each side of the footrest
- Foot section release handles operate with one hand
- Stainless-steel debris pan slides out from the end of the footrest
- Chair arms rotate completely out of the way for chair accessibility
- Powered multi-positioning chair adjusts the seat, height, back, and tilt for a more ergonomic working position
- Back support angles from 0° to 80°
- Seat tilt angles from 0° to 30°
- Low 21" seat height and wide, flat surface enable easier patient transfer on and off the chair
- Premium comfort system uses unique, contoured upholstery to help cradle the patient and reduce movement during procedures
- Seamless upholstery sanitizes easier for efficient patient turnover (includes a plastic protective foot cover)

### Specifications:

- Length: 75" (84" with footrest extended)
- Seat Section Width: 29" (74 cm)
- Foot Section Width: 26" (66 cm)
- Seat Height: 19" (minimum) - 40" (maximum)
- Back Angle Adjustment: 0 – 80 Degrees
- Load Capacity: 450 lbs (204 kg)
- Upholstery: Premium Seamless or UltraFree (looks and feels like leather)
- Duplex, Hospital-grade outlets (on each side)
- Four customizable positions

*\*\*Brochure and images are a general overview of the product. Actual product options may vary.*

**VENTUREMEDICAL.COM**