



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Midmark 419 ENT Chair Refurbished



- Modern Design
- Narrow Headrest
- Slim Backrest
- Foot Control
- Electrically Operated
- +/- Height, +/- Backrest
- +/- Foot, +/- Tilt
- 22.5" Low Height
- Adjustable, Articulating Headrest
- Auto-Return
- Armrests
- Completely Refurbished

The Midmark 419 Power Exam Chair is ideal for Dermatologists, Plastic Surgeons, Oral or Maxillofacial Surgeons and all healthcare professionals requiring maximum access to a patient's head and neck. It features a narrow neck area and a unique, slim backrest, for a large work surface and greater legroom for improved patient access and comfort.

The Midmark 419 Power Exam Chair powers all the way down to 22 3/8" so patients can easily get on and off the chair, without the need to be lifted. Once the patient is on the chair, it powers up to as high as 40 3/8", so physicians won't have to bend or stoop to perform their procedure.

The Midmark 419 comes completely refurbished with a comprehensive warranty and new upholstery.

Specifications:

- Minimum Height: 22.5" (55.9cm)
- Maximum Height: 40.5" (78.7cm)
- Trendelenburg to 30°
- Width: 23" (61cm)
- Length (extended): 65" (165.1cm)
- Max Patient Weight: 325 lbs
- Electrical: Domestic: 115v, 11AMP