

## Midmark 222 Power Exam Table Refurbished



Stock Photo – Options May Vary

The Midmark 222 Power exam table is one of the most patient accessible exam tables in the industry. The table lowers to a low height of 18" and lays flat for easy and safe patient entry and exit. The table features pull out stumps, power height and backrest adjustment, extendable footrest with pad, hospital grade power outlet and a paper roll holder.

The Midmark 222 Exam Table also includes a pull-out drain pan at the end of the seat section and pull-out heel stumps.

The Midmark 222 comes completely refurbished with new upholstery (your choice of colors), foot control, table paper roll holder, and a comprehensive warranty.

- Maximum Patient Accessibility
- Power height and Back
- Foot Control
- Low 18" Minimum Height
- Weight Capacity (400 lbs)
- Pull Out Stirrups
- Extendable Leg Section with Pad
- Paper Roll Holder
- Table Paper Holder
- Table Paper Strap
- Completely Refurbished

### Specifications:

- Weight Capacity: 400 lbs
- Tabletop Height: 18"- 37"
- Pelvic Tilt: 7°
- Dimensions:
  - Upholstered Top: 28.5" W x 59.5" L
  - With Footrest Extended: 74" L
- Paper Rolls: up to 21" x 3.5"
- Power: 115 VACS, 5 A Max

*\*\*Brochure and images are a general overview of the product.  
Actual product options may vary.*