

## Alaris Medley 8015 Pump w/ 8100 Infusion Module - Refurbished



Stock Photo – Options May Vary

- Battery Operation
- Patient ID Entry
- Guardrails Safety Software
- Concentration Limits Data Set
- Anesthesia Mode
- Audible Alarm
- Dose Checking
- Electric Memory
- Completely Refurbished

The Alaris Medley medication safety system is a modular infusion and monitoring system intended for use in today's growing professional system. The Point of Care Unit (8000 Series) consists of the 8015 core, Guardrails Safety Software, and up to 4 detachable modules which provide infusion or monitoring capabilities.

The Alaris 8015 system can be used for adult, pediatric, and neonatal care. The Guardrail Safety Software brings a new level of medication error prevention to the point of patient care.

The Alaris Medley 8015 comes completely refurbished with a comprehensive warranty.

### Specifications:

- Overall Dims (8000PC):
  - 6.9"W x 8.8"H x 9"D
- Overall Dims (8100):
  - 3.3"W x 8.9"H x 5.5"D
- Battery: 8 hrs w/ 1 Pump at 25mL/h
- Data Port: RS-232 w/ RJ45
- Weight (8000PC): 7.2 lbs
- Weight (8100): 2.5 lbs
- Power Requirements:
  - 100-240V, 150VA Max
  - 50/60Hz

*\*\*Brochure and images are a general overview of the product. Actual product options may vary.*