

Drager Fabius GS Anesthesia Workstation - Refurbished



Stock Photo – Options May Vary

The Fabius GS combines advanced ventilation technology with reliable and proven breathing system components, opening new frontiers in ventilation performance while maintaining well-established functionality.

The Fabius GS is suitable for any patient, pediatric to adult. Pressure Support mode facilitates spontaneous breathing by removing the work of breathing due to circuit resistance, improving comfort levels and enhancing quality of patient care.

The modern technology and conveniences of the Drager Fabius GS, coupled with the upgradeable features and options makes this unit a smart choice for any anesthesia setting. Optional heated breathy circuit (Cosy), SIMV and Pressure Support available for an additional cost.

The Drager Fabius GS comes completely refurbished with a comprehensive warranty.

- Modern Technology
- Extremely Safe for Patients
- Integrated Absorber
- O₂ and Tidal Volume Monitoring
- Integrated Breathing Circuit
- Electronic Ventilator
- Electronic PEEP
- Pressure and Volume Control
- AGSS Scavenging System w/ Sample Gas Return
- Low Flow Capable
- 2 Vaporizer Manifold
- Storage Drawers
- Compact Unit for places where space is a premium
- Built-in back-up battery to handle sudden power-down situations
- Completely Refurbished

Specifications:

- Height (unit only): 51.2"
- Width (unit only): 35.2"
- Depth (unit only): 32.3"
- Weight (unit only): 224 Lbs
- Power Requirement: 120v, 2.3A

***Brochure and images are a general overview of the product.
Actual product options may vary.*