

# Physio-Control LifePak 20e Defibrillator Biphasic | Pacing | AED - Refurbished



Stock Photo – Options May Vary

The Lifepak 20e defibrillator/monitor is the ideal crash cart defibrillator, putting early, effective defibrillation into the hands of BLS responders while transitioning easily to full therapy and monitoring capabilities for the ALS trained code team.

The Lifepak 20e will come with Biphasic waveform technology, non-invasive pacing and AED mode.

Included with the unit: Multi-Function Cable, ECG Cable and Lead wires, Set of Adult Therapy Pads, Printer Paper and New Battery.

Optional accessories and functions include: SpO<sub>2</sub>, Carrying Case and External Paddles.

The Lifepak 20e defibrillator comes completely refurbished with a comprehensive warranty.

- Biphasic Waveform Technology
- Pacing Function
- AED Function
- Synchronized Cardioversion
- 3 Lead EKG
- Onboard Thermal Printer
- Hands Free Multi-Function Cable
- Built-In Battery and Charger
- Small, Lightweight Design
- Carry Handle
- Many Optional Functions Available
- Completely Refurbished

## Specifications:

- Max Joules: 360j
- Joule Range: 2-360
- Charge Time: Less than 7 seconds
- Pacing Type: Demand or Non-Demand
- Display Size (Screen): 4.53”W x 3.4”H
- Battery: Rechargeable, Lithium-Ion
- Recorder Paper: 50mm, Thermal
- Height (unit only): 21.3 cm (8.4”)
- Width (unit only): 26.2 cm (10.3”)
- Depth (unit only): 26.2 cm (10.3”)
- Weight (unit only): ~5.58 kg (~12.3 lbs)
- Power Requirement: 120v

*\*\*Brochure and images are a general overview of the product.  
Actual product options may vary.*