

Zoll R-Series Defibrillator - Refurbished



Stock Photo – Options May Vary

- Biphasic Waveform Technology
- Pacing Function
- Synchronized Cardioversion
- 3 Lead EKG
- Onboard Thermal Printer
- Hands Free Multi-Function Cable
- Built-In Battery Charger
- One Step Cable Manager
- See-Thru CPR
- Many Optional Functions Available
- Monitor Displays Selected and Delivered Energy
- Completely Refurbished

The Zoll R-Series defibrillator is designed to help you deliver optimal care to adult and pediatric patients. This unit features biphasic waveform, non-invasive pacing and see-thru CPR. Operating the system is simple and user-friendly with straightforward controls and easy to use menus.

Included with the unit: Multi-Function Cable, ECG Cable and Leadwires, Set of Adult Therapy Pads, Printer Paper and New Battery.

Optional accessories and functions include: 5 Lead ECG, SpO₂, NIBP, ETCO₂ (Mainstream or Sidestream), and External Paddles.

The Zoll R-Series comes completely refurbished with a comprehensive warranty.

Specifications:

- Max Joules: 200j
- Charge Time: Less than 7 seconds
- Display Size (Screen): 6.5"
- Battery: Rechargeable, Lithium Ion
- Recorder Paper: 90mm, Thermal
- Height (unit only): 8.2"
- Width (unit only): 10.5"
- Depth (unit only): 12.5"
- Weight (unit only): 13.6 lbs w/o Paddles
- Power Requirement: 100-120V AC

***Brochure and images are a general overview of the product. Actual product options may vary.*