



www.venturemedical.com

Venture Medical ReQuip, Inc.
6008 Bonacker Drive
Tampa, FL 33610
800-627-3215 | 813-740-8640
Fax: 800-466-1251
sales@venturemedical.com

Stryker 1115 Prime Big Wheel Stretcher -Refurbished



Stock Photo – Options May Vary

The Stryker 1115 Prime Series Stretchers are extraordinarily mobile and easy to use. Big Wheel steering helps make transport safe and efficient while reducing the physical strain on the clinician.

Features include: Dual-side mounted foot controls with uni-lower pedal, dual end brake/steer controls, glideaway siderails with integrated push/pull handles, pneumatic backrest.

The Stryker 1115 Prime stretcher comes completely refurbished with a comprehensive warranty and a new pad.

- Glideaway Side Rails
- Central Brake and Steer
- $\pm$ Height, $\pm$ Trendelenburg
- Pneumatic Backrest
- O2 Tank Storage
- 700 lbs Weight Capacity
- IV Pole (optional)
- Includes New 4" Pad
- Completely Refurbished

Specifications:

- Overall Length: 85"
- Overall Width: 34"
- Patient Surface Width: 26"
- Height Range: 20.75"– 34"
- Weight Capacity: 700 lbs
- Fowler Angle: 0° - 90°
- Trend/Rev: +16° to -17°
- Wheel Diameter: 8"
- Pad Thickness: 4"

***Brochure and images are a general overview of the product.
Actual product options may vary.*

VENTUREMEDICAL.COM

© 2007-2026 Venture Medical ReQuip, Inc.