

Stryker 735 Transport Stretcher - Refurbished



Stock Photo – Options May Vary

- Gas Spring Assist Backrest
- Fold Down Side Rails
- Central Brake and Steer
- ± Height
- ± Trendelenburg
- Built-in O2 Tank Holder
- 500 lbs Weight Capacity
- Includes New 4" Pad
- Completely Refurbished

The Stryker 735 Stretcher has a simple design that provides safe, efficient transport within the hospital, or surgery center. The 5th wheel provides easier steering and cornering.

Additional features include: Gas-spring assisted backrest (fowler), trendelenburg, reverse trendelenburg, and height adjustment. The 500 lbs weight capacity can handle large patients with ease.

The Stryker 735 stretcher comes completely refurbished with a comprehensive warranty and a new pad.

Specifications:

- Overall Length: 83"
- Overall Width: 30"
- Patient Surface Width: 26"
- Height Range: 21.5" – 36"
- Weight Capacity: 500 lbs
- Trend/Rev Trend: ± 18°
- Side Rails: 13" x 55"
- Wheel Diameter: 8"
- Pad Thickness: 4"

***Brochure and images are a general overview of the product.
Actual product options may vary.*